

**OPTIONS**


**Lock and Key Options**

Option Code	Description
30	Standard locks on doors/drawers supplied and factory installed *
31	Provisions for standard locks on doors/drawers (holes & tabs) *
37	Hasp Lock (padlocks not included)
37A	Hasp lock with stainless steel plate (pad lock not included)
38	Customer specified locks on doors/drawers supplied and factory installed **
39	Customer specified locks on doors/drawers supplied and not factory installed **
BS	Best locks (housing only), core not provided
LK4	Standard 5 pin tumbler lock with removable core, three level keying

Notes: \*Standard lock is a disc tumbler lock. Each cabinet is keyed differently per unit unless otherwise specified.


\*\*Provide complete specification and source for all special options ordered.

**All Locks Keyed Different Per Unit**




Locks are keyed different per unit from other units.  
Note: This is the standard key schedule with disc tumbler lock supplied for option code 30 unless otherwise specified.

**All Locks Keyed Alike**



All locks are keyed alike, each lock can be opened with a single key.

**Locks Keyed Alike in Sets**




Locks are keyed alike in sets. Only one key will open one set of locks within the set. Up to 229 key changes for disc are available.




Locks are keyed different from other locks. Up to 229 key changes for disc are available.

**Master Keyed Locks**



All locks are keyed different from all other locks. Each set within a group can be opened with a master key. Up to 229 key changes for disc are available.

**Pin Locks Keyed Different, Grand Master & Master Keyed**



All locks are keyed different from all other locks. Each set within a group can be opened with a master key. All groups can be opened with a grand master key. 400 key changes and up to 10 master keys.

**OPTIONS**

**Lock & Key Options**



- Option 30 (Standard)\*  
5-Disc Master,  
removable core
- Option IK4  
5-Pin Tumbler, removable  
core, 3 level keying




- Option 37  
Hasp Lock  
(padlock not included)



- Option 37A  
Hasp Lock with stainless  
steel plate  
(padlock not included)


- Option 31 - Provisions for standard locks on doors/drawers (holes and tabs)\*
  - Option 38 - Customer specified locks on doors/drawers supplied and factory installed\*\*
  - Option 39 - Customer specified locks on doors/drawers supplied and not factory installed \*\*
  - Option BS - Best locks (housing only), core not provided
- \*Standard lock is a disc tumbler lock. Each cabinet is keyed differently per unit unless otherwise specified.  
\*\*Provide complete specification and source for all special options ordered.

**All Locks Keyed Different Per Unit**




Locks are keyed different per unit from other units.  
Note: This is the standard key schedule with disc tumbler lock supplied for option code 30 unless otherwise specified.

**Master Keyed Locks**




All locks are keyed different from all other locks. Each set within a group can be opened with a master key. Up to 229 key changes for disc are available.

**All Locks Keyed Alike**




All locks are keyed alike, each lock can be opened with a single key.

**Pin Locks Keyed Different, Grand Master & Master Keyed**




All locks are keyed different from all other locks. Each set within a group can be opened with a master key. All groups can be opened with a grand master key. 400 key changes and up to 10 master keys.

**Locks Keyed Alike in Sets**



Locks are keyed alike in sets. Only one key will open one set of locks within the set. Up to 229 key changes for disc are available.



Locks are keyed different from other locks. Up to 229 key changes for disc are available.

All dimensions and sizes shown are nominal. Specifications and details are based on product information at the time of printing and may change at any time without notice. Mott Manufacturing reserves the right to change dimensions, specifications and manufacturing details at any time without notice.