

the right formula

Amherst College - Geology & Natural History Museum

OVERVIEW

- **Project:** 60,000 sq. ft. earth sciences research and education facility
- **Architect:** Payette
- **GC:** Daniel O'Connell's Sons





In this new space for Amherst College - Geology & Natural History Museum

Lexington Series Wood Veneer Casework

Epoxy Resin Tops

Harvard Series Fume Hoods

The 60,000 sq. ft. earth sciences building and museum includes teaching labs and specialized research spaces. This project features Lexington Series wood casework with quartered cut white oak veneer, with a dark stain and black epoxy work surfaces. Custom-designed tall-storage were built to accommodate the school's unique rock storage trays. Also featured are several Harvard Series fume hoods.