

the right formula

Brown University - Biomed Labs

OVERVIEW

- **Project:** 5-story university and research center
- **Architect:** Tsoi/Kobus & Associates Inc.
- **GC:** Shawmut Design & Construction





In this new space for Brown University - Biomed Labs

Wakefield Series Design Overlay Casework

Epoxy Resin and Plastic Laminate Tops

Epoxy Resin Sinks and Pegboards

Pro Series Fume Hoods

Brown University was designed with five floors of Wakefield Series casework constructed of painted steel frames and complemented with vertically matched maple door and drawer fronts. Sitting height writing stations with plastic laminate countertops and shelves with tack boards, were included at each bench. Black epoxy counter tops, sinks and pegboards included with WaterSaver faucets and fixtures. Custom cabinets were built with extra wide drawers ideally used for flat storage of maps, posters or charts. Finally, Pro Series fume hoods which contain a chain and sprocket sash system were also supplied on this project.