

the right formula for Coca-Cola Bottling Co. - Testing Lab

OVERVIEW

- **Project:** Food & beverage industrial laboratory
- **Architect:** CMC Associates



new england lab[®]
lab-tested furniture systems



In this new lab for Coca-Cola Bottling Co. - Testing Lab

[Manchester Series Painted Steel Casework](#)

[Wall Cabinets with Custom Green Colored Glass](#)

[Black Epoxy Resin Tops](#)

[Pro Series Fume Hoods](#)

This food & beverage industrial lab features Manchester Series painted steel casework, black epoxy resin work surfaces and Pro Series bench fume hoods. Wall cabinets were fit with a custom green colored glass and many of them were suspended from the ceiling in a pass-through configuration. Tilt out waste bins were also implemented in this testing lab.