


FUME HOODS - HARVARD SERIES

Bench Hood - Frameless Vertical Sash




- Features:**
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -FRP Liner
 -54-1/4" High

| Width | Item Number |
|-------|-------------|
| 36" | 7111070 |
| 48" | 7211070 |
| 60" | 7311070 |
| 72" | 7411070 |
| 96" | 7511070 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |

Bench Hood - Combination Sash




- Features:**
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -FRP Liner
 -54-1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 7211050 |
| 60" | 7311050 |
| 72" | 7411050 |
| 96" | 7510050 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-----------------------|------|-----------------------|------|---------------------|------|---------------------|------|
| | | 1/2 SLIDING DOOR OPEN | | 1/2 SLIDING DOOR OPEN | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 48" | 10" | 315 | 0.05 | 394 | 0.05 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 414 | 0.05 | 517 | 0.05 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 512 | 0.05 | 640 | 0.10 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 709 | 0.05 | 886 | 0.10 | 913 | 0.07 | 1142 | 0.11 |

FUME HOODS - HARVARD SERIES

Bench Hood - Horizontal Sash




Features:
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Framed Horizontal Sliding Sash
 -FRP Liner
 -54-1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 7211020 |
| 60" | 7311020 |
| 72" | 7411020 |
| 96" | 7511020 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|-----------------------|------|-----------------------|------|
| | | 1/2 SLIDING DOOR OPEN | | 1/2 SLIDING DOOR OPEN | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 48" | 10" | 315 | 0.05 | 394 | 0.05 |
| 60" | 12" | 414 | 0.05 | 517 | 0.05 |
| 72" | 12" | 512 | 0.05 | 640 | 0.10 |
| 96" | 2@10" | 709 | 0.05 | 886 | 0.10 |

Bench Hood - Framed Vertical Sash




Features:
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -FRP Liner
 -54-1/4" High

| Width | Item Number |
|-------|-------------|
| 36" | 7111060 |
| 48" | 7211060 |
| 60" | 7311060 |
| 72" | 7411060 |
| 96" | 7511060 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |

FUME HOODS - HARVARD SERIES

Bench Hood - Fixed Center Post




- Features:**
- 34 5/8" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Side Access Panels
 - FRP Liner
 - 54 1/4" High

| Width | Item Number |
|-------|-------------|
| 72" | 7411090 |
| 96" | 7511090 |
| 120" | 7611090 |
| 144" | 7H11090 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 Max SASH OPENING | | 28-1/2 Max SASH OPENING | | 18 Max SASH OPENING | | 18 Max SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |
| 120" | 2@12" | 1811 | 0.15 | 2264 | 0.20 | 1167 | 0.06 | 1458 | 0.09 |
| 144" | 2@12" | 2205 | 0.20 | 2755 | 0.30 | 1420 | 0.08 | 1775 | 0.13 |

Bench Hood - Disappearing Center Post




- Features:**
- 34 5/8" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Interior Side Access Panels
 - FRP Liner
 - 54 1/4" High

| Width | Item Number |
|-------|-------------|
| 72" | 7411080 |
| 96" | 7511080 |
| 120" | 7611080 |
| 144" | 7411080 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 Max SASH OPENING | | 28-1/2 Max SASH OPENING | | 18 Max SASH OPENING | | 18 Max SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |
| 120" | 2@12" | 1811 | 0.15 | 2264 | 0.20 | 1167 | 0.06 | 1458 | 0.09 |
| 144" | 2@12" | 2205 | 0.20 | 2755 | 0.30 | 1420 | 0.08 | 1775 | 0.13 |

FUME HOODS - HARVARD SERIES

Floor-Mounted Hood - Frameless Vertical Sash




Features:
 -38" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Side Access Panels
 -FRP Liner
 -54 1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 7211071 |
| 60" | 7311071 |
| 72" | 7411071 |
| 96" | 7511071 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 |

Floor-Mounted Hood - Combination Sash




Features:
 -38" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Side Access Panels
 -FRP Liner
 -54 1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 7211051 |
| 60" | 7311051 |
| 72" | 7411051 |
| 96" | 7511051 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 |

FUME HOODS - HARVARD SERIES

Floor-Mounted Hood - Horizontal Sash




- Features:**
- 38" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Interior Side Access Panels
 - FRP Liner
 - 90" High

| Width | Item Number |
|-------|-------------|
| 48" | 7211021 |
| 60" | 7311021 |
| 72" | 7411021 |
| 96" | 7511021 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|------------------------|------|------------------------|------|
| | | 1/2 SLIDING DOORS OPEN | | 1/2 SLIDING DOORS OPEN | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 48" | 10" | 680 | 0.20 | 849 | 0.25 |
| 60" | 12" | 890 | 0.15 | 1113 | 0.20 |
| 72" | 12" | 1102 | 0.20 | 1378 | 0.35 |
| 96" | 2@10" | 1525 | 0.20 | 1907 | 0.35 |

Floor-Mounted Hood - Vertical Sash




- Features:**
- 38" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Interior Side Access Panels
 - FRP Liner
 - 90" High

| Width | Item Number |
|-------|-------------|
| 48" | 7211061 |
| 60" | 7311061 |
| 72" | 7411061 |
| 96" | 7511061 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 |

FUME HOODS - HARVARD SERIES

Floor-Mounted Hood - Double Depth - Combination Sash




- Features:**
- 69 1/4" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Interior Side Access Panels
 - FRP Liner
 - 90" High

| Width | Item Number |
|-------|-------------|
| 72" | 7411057 |
| 96" | 7511057 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 |

Floor-Mounted Hood - Double Sided, Double Depth - Horizontal Sash




- Features:**
- 69 1/4" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Interior Side Access Panels
 - FRP Liner
 - 90" High

| Width | Item Number |
|-------|-------------|
| 72" | 7411026 |
| 96" | 7511026 |
| 120" | 7611026 |
| 144" | 7H11026 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|------------------------|------|------------------------|------|
| | | 1/2 SLIDING DOORS OPEN | | 1/2 SLIDING DOORS OPEN | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 72" | 12" | 1102 | 0.20 | 1378 | 0.35 |
| 96" | 2@10" | 1525 | 0.20 | 1907 | 0.35 |
| 120" | 2@12" | 1950 | 0.15 | 2436 | 0.25 |
| 144" | 2@12" | 2373 | 0.25 | 2965 | 0.35 |

FUME HOODS HARVARD SERIES

Bench Hood - ADA Height - Combination Sash




Features:
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Side Access Panels
 -FRP Liner
 -66 1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 72B1010 |
| 60" | 73B1010 |
| 72" | 74B1010 |
| 96" | 75B1010 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |

Radioisotope Bench Hood - Combination Sash




Features:
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -316 S/S Radiused Corners
 -54 1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 7218050 |
| 60" | 7318050 |
| 72" | 7418050 |
| 96" | 7518050 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-----------------------|------|-----------------------|------|---------------------|------|---------------------|------|
| | | 1/2 SLIDING DOOR OPEN | | 1/2 SLIDING DOOR OPEN | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 48" | 10" | 315 | 0.05 | 394 | 0.05 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 414 | 0.05 | 517 | 0.05 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 512 | 0.05 | 640 | 0.10 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 709 | 0.05 | 886 | 0.10 | 913 | 0.07 | 1142 | 0.11 |

FUME HOODS HARVARD SERIES

Radioisotope Bench Hood - Framed Vertical Sash




- Features:**
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -316 S/S Radiused Corners
 -54 1/4" High

| Width | Item Number |
|-------|-------------|
| 36" | 7118060 |
| 48" | 7218060 |
| 60" | 7318060 |
| 72" | 7418060 |
| 96" | 7518060 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |

Perchloric Acid Bench Hood - Combination Sash




- Features:**
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -316 S/S Radiused Corners
 -54 1/4" High

| Width | Item Number |
|-------|-------------|
| 48" | 7217050 |
| 60" | 7317050 |
| 72" | 7417050 |
| 96" | 7517050 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-----------------------|------|-----------------------|------|---------------------|------|---------------------|------|
| | | 1/2 SLIDING DOOR OPEN | | 1/2 SLIDING DOOR OPEN | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 48" | 10" | 315 | 0.05 | 394 | 0.05 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 414 | 0.05 | 517 | 0.05 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 512 | 0.05 | 640 | 0.10 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 709 | 0.05 | 886 | 0.10 | 913 | 0.07 | 1142 | 0.11 |

FUME HOODS - HARVARD SERIES

Perchloric Acid Bench Hood - Vertical Framed Sash




- Features:**
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -316 S/S Radiused Corners
 -54 1/4" High

| Width | Item Number |
|-------|-------------|
| 36" | 7117060 |
| 48" | 7217060 |
| 60" | 7317060 |
| 72" | 7417060 |
| 96" | 7517060 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 28-1/2 MAX SASH OPENING | | 28-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 434 | 0.05 | 543 | 0.10 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 630 | 0.15 | 789 | 0.20 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 830 | 0.15 | 1035 | 0.20 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1025 | 0.20 | 1280 | 0.25 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1418 | 0.20 | 1772 | 0.25 | 913 | 0.07 | 1142 | 0.11 |

Extra High Bench Hood - Frameless Vertical Sash




- Features:**
 -66 1/4" High
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -FRP Liner

| Width | Item Number |
|-------|-------------|
| 36" | 7111075 |
| 48" | 7211075 |
| 60" | 7311075 |
| 72" | 7411075 |
| 96" | 7511075 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|------|---------------------|------|---------------------|------|
| | | 32-1/2 MAX SASH OPENING | | 32-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 479 | 0.08 | 598 | 0.12 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 695 | 0.17 | 869 | 0.25 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 912 | 0.14 | 1140 | 0.21 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1129 | 0.21 | 1411 | 0.32 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1562 | 0.21 | 1952 | 0.32 | 913 | 0.07 | 1141 | 0.11 |

FUME HOODS - HARVARD SERIES

Extra High Bench Hood - Combination Sash




- Features:**
 -66 1/4" High
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -FRP Liner

| Width | Item Number |
|-------|-------------|
| 48" | 7211055 |
| 60" | 7311055 |
| 72" | 7411055 |
| 96" | 7511055 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|---------------------|------|---------------------|------|---------------------|------|---------------------|------|
| | | 50 MAX SASH OPENING | | 50 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 239 | 0.04 | 299 | 0.06 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 347 | 0.08 | 434 | 0.12 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 456 | 0.07 | 570 | 0.11 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 564 | 0.11 | 705 | 0.16 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 781 | 0.11 | 976 | 0.16 | 913 | 0.07 | 1141 | 0.11 |

Extra High Bench Hood - Horizontal Sash




- Features:**
 -66 1/4" High
 -34 5/8" Deep
 -Hinged Light Fixture
 -Vision Panel
 -Removable Interior Side Access Panels
 -FRP Liner

| Width | Item Number |
|-------|-------------|
| 48" | 7211025 |
| 60" | 7311025 |
| 72" | 7411025 |
| 96" | 7511025 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | |
|--------------------|----------|---------------------|------|---------------------|------|
| | | 50 MAX SASH OPENING | | 50 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP |
| 36" | 10" | 239 | 0.04 | 299 | 0.06 |
| 48" | 10" | 347 | 0.08 | 434 | 0.12 |
| 60" | 12" | 456 | 0.07 | 570 | 0.11 |
| 72" | 12" | 564 | 0.11 | 705 | 0.16 |
| 96" | 2@10" | 781 | 0.11 | 976 | 0.16 |

FUME HOODS - HARVARD SERIES

Extra High Bench Hood - Framed Vertical Sash



- Features:**
- 66 1/4" High
 - 34 5/8" Deep
 - Hinged Light Fixture
 - Vision Panel
 - Removable Interior Side Access Panels
 - FRP Liner

| Width | Item Number |
|-------|-------------|
| 36" | 7111065 |
| 48" | 7211065 |
| 60" | 7311065 |
| 72" | 7411065 |
| 96" | 7511065 |

| EXHAUST PARAMETERS | | 80 FPM | | 100 FPM | | 80 FPM | | 100 FPM | |
|--------------------|----------|-------------------------|------|-------------------------|-------|---------------------|------|---------------------|------|
| | | 32-1/2 MAX SASH OPENING | | 32-1/2 MAX SASH OPENING | | 18 MAX SASH OPENING | | 18 MAX SASH OPENING | |
| HOOD SIZE | DUCT DIA | CFM | SP | CFM | SP | CFM | SP | CFM | SP |
| 36" | 10" | 479 | 0.08 | 598 | 0.012 | 280 | 0.03 | 350 | 0.04 |
| 48" | 10" | 695 | 0.17 | 869 | 0.25 | 407 | 0.06 | 508 | 0.09 |
| 60" | 12" | 912 | 0.14 | 1140 | 0.21 | 533 | 0.05 | 667 | 0.07 |
| 72" | 12" | 1129 | 0.21 | 1411 | 0.32 | 660 | 0.07 | 825 | 0.11 |
| 96" | 2@10" | 1562 | 0.21 | 1952 | 0.32 | 913 | 0.07 | 1141 | 0.11 |