

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Note: Add option code 15 to order a finished back on these units.

Open Units



| Width | Item Number |
|--------|-------------|
| 15" | 1950000 |
| 18" | 1150000 |
| 24" | 1250000 |
| 30" | 1350000 |
| 36" | 1450000 |
| 42" | 1850000 |
| 48" | 1550000 |
| 513mm | 1650000 |
| 1000mm | 1750000 |

Full Door Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1950011 | 1950211 | 1 |
| 18" | 1150011 | 1150211 | 1 |
| 24" | 1250011 | 1250211 | 1 |
| 30" | 1350011 | - | 2 |
| 36" | 1450011 | - | 2 |
| 42" | 1850011 | - | 2 |
| 48" | 1550011 | - | 2 |
| 513mm | 1650011 | 1650211 | 1 |
| 1000mm | 1750011 | - | 2 |

- Cabinets receive one shelf.

Door/Drawer Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1950022 | 1950222 | 1 |
| 18" | 1150022 | 1150222 | 1 |
| 24" | 1250022 | 1250222 | 1 |
| 30" | 1350022 | - | 2 |
| 36" | 1450022 | - | 2 |
| 42" | 1850022 | - | 2 |
| 48" | 1550022 | - | 2 |
| 513mm | 1650022 | 1650222 | 1 |
| 1000mm | 1750022 | - | 2 |

- Cabinets up to 24" wide receive one shelf.
- Cabinets 25" and wider do not have a shelf.



| Width | Right Hand | Left Hand |
|--------|------------|-----------|
| 30" | 1350041 | 1350014 |
| 36" | 1450041 | 1450014 |
| 42" | 1850041 | 1850014 |
| 48" | 1550041 | 1550014 |
| 1000mm | 1750041 | 1750014 |

- Left hand Door/Drawer unit shown. Right hand Door/Drawer units are reverse of what is shown.
- Cabinet receives one shelf.

Sliding Door Units



| Width | Item Number |
|--------|-------------|
| 30" | 1350711 |
| 36" | 1450711 |
| 42" | 1850711 |
| 48" | 1550711 |
| 1000mm | 1750711 |

- Cabinets receive one shelf.

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Note: Add option code 15 to order a finished back on these units.

Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950055 |
| 18" | 1150055 |
| 24" | 1250055 |
| 30" | 1350055 |
| 36" | 1450055 |
| 42" | 1850055 |
| 48" | 1550055 |
| 513mm | 1650055 |
| 1000mm | 1750055 |



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950066 |
| 18" | 1150066 |
| 24" | 1250066 |
| 30" | 1350066 |
| 36" | 1450066 |
| 42" | 1850066 |
| 48" | 1550066 |
| 513mm | 1650066 |
| 1000mm | 1750066 |



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950044 |
| 18" | 1150044 |
| 24" | 1250044 |
| 30" | 1350044 |
| 36" | 1450044 |
| 42" | 1850044 |
| 48" | 1550044 |
| 513mm | 1650044 |
| 1000mm | 1750044 |



| Width | Item Number |
|--------|-------------|
| 30" | 1350097 |
| 36" | 1450097 |
| 42" | 1850097 |
| 48" | 1550097 |
| 1000mm | 1750097 |



| Width | Item Number |
|--------|-------------|
| 30" | 1350098 |
| 36" | 1450098 |
| 42" | 1850098 |
| 48" | 1550098 |
| 1000mm | 1750098 |

File Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950122 |
| 18" | 1150122 |
| 24" | 1250122 |
| 30" | 1350122 |
| 36" | 1450122 |
| 42" | 1850122 |
| 48" | 1550122 |
| 513mm | 1650122 |
| 1000mm | 1750122 |

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Note: Add option code 15 to order a finished back on these units.

File Drawer Units



up to 24" wide



30" and wider



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950177 |
| 18" | 1150177 |
| 24" | 1250177 |
| 30" | 1350177 |
| 36" | 1450177 |
| 42" | 1850177 |
| 48" | 1550177 |
| 513mm | 1650177 |
| 1000mm | 1750177 |

| Width | Item Number |
|--------|-------------|
| 15" | 1950188 |
| 18" | 1150188 |
| 24" | 1250188 |
| 30" | 1350188 |
| 36" | 1450188 |
| 42" | 1850188 |
| 48" | 1550188 |
| 513mm | 1650188 |
| 1000mm | 1750188 |

Sink Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1953011 | 1953211 | 1 |
| 18" | 1153011 | 1153211 | 1 |
| 24" | 1253011 | 1253211 | 1 |
| 30" | 1353011 | - | 2 |
| 36" | 1453011 | - | 2 |
| 42" | 1853011 | - | 2 |
| 48" | 1553011 | - | 2 |
| 513mm | 1653011 | 1653211 | 1 |
| 1000mm | 1753011 | - | 2 |

- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet.
- Add option code 90 to order Sink Hanger Supports.

Acid Storage Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1951011 | 1951211 | 1 |
| 18" | 1151011 | 1151211 | 1 |
| 24" | 1251011 | 1251211 | 1 |
| 30" | 1351011 | - | 2 |
| 36" | 1451011 | - | 2 |
| 42" | 1851011 | - | 2 |
| 48" | 1551011 | - | 2 |
| 513mm | 1651011 | 1651211 | 1 |
| 1000mm | 1751011 | - | 2 |

- Cabinet receives a molded liner and one chemically resistant shelf.
- Add option code VK to order a Corrugated Polyolefin Vent Kit for these units.
- Note: Acid Cabinets must be vented. Mott recommends venting acid cabinets with a Polyolefin Vent Kit (PVSA for FRP lined cabinets).

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Note: Add option code 15 to order a finished back on these units.

Solvent Storage Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1952211 | 1952311 | 1 |
| 18" | 1152211 | 1152311 | 1 |
| 24" | 1252211 | 1252311 | 1 |
| 30" | 1352211 | - | 2 |
| 36" | 1452211 | - | 2 |
| 42" | 1852211 | - | 2 |
| 48" | 1552211 | - | 2 |
| 513mm | 1652211 | 1652311 | 1 |
| 1000mm | 1752211 | - | 2 |

- Solvent units are not insulated.
- Add option code 83 for Self Closing/Latching doors.
- Add option code "FL" for Self Closing/Latching doors with a fusible link.
- Solvent Units are designed to meet OSHA and NFPA 30 standards.
- See page E17 for Solvent cabinet details.

SUSPENDED CABINETS

18-3/4" HIGH X 22" DEEP

Note: Add option code 15 to order a finished back on these units.

Open Units



| Width | Item Number |
|--------|-------------|
| 15" | 1960000 |
| 18" | 1160000 |
| 24" | 1260000 |
| 30" | 1360000 |
| 36" | 1460000 |
| 42" | 1860000 |
| 48" | 1560000 |
| 513mm | 1660000 |
| 1000mm | 1760000 |

Full Door Unit



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1960011 | 1960211 | 1 |
| 18" | 1160011 | 1160211 | 1 |
| 24" | 1260011 | 1260211 | 1 |
| 30" | 1360011 | - | 2 |
| 36" | 1460011 | - | 2 |
| 42" | 1860011 | - | 2 |
| 48" | 1560011 | - | 2 |
| 513mm | 1660011 | 1660211 | 1 |
| 1000mm | 1760011 | - | 2 |

Door/Drawer Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1960022 | 1960222 | 1 |
| 18" | 1160022 | 1160222 | 1 |
| 24" | 1260022 | 1260222 | 1 |
| 30" | 1360022 | - | 2 |
| 36" | 1460022 | - | 2 |
| 42" | 1860022 | - | 2 |
| 48" | 1560022 | - | 2 |
| 513mm | 1660022 | 1660222 | 1 |
| 1000mm | 1760022 | - | 2 |



| Width | Item Number |
|--------|-------------|
| 30" | 1360041 |
| 36" | 1460041 |
| 42" | 1860041 |
| 48" | 1560041 |
| 1000mm | 1760041 |

- Right hand Door/Drawer Unit shown. Left hand reversed.

SUSPENDED CABINETS

18-3/4" HIGH X 22" DEEP

Note: Add option code 15 to order a finished back on these units.

Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1960044 |
| 18" | 1160044 |
| 24" | 1260044 |
| 30" | 1360044 |
| 36" | 1460044 |
| 42" | 1860044 |
| 48" | 1560044 |
| 513mm | 1660044 |
| 1000mm | 1760044 |

File Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1960177 |
| 18" | 1160177 |
| 24" | 1260177 |
| 30" | 1360177 |
| 36" | 1460177 |
| 42" | 1860177 |
| 48" | 1560177 |
| 513mm | 1660177 |
| 1000mm | 1760177 |

Sink Units



up to 24" wide




30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1963011 | 1963211 | 1 |
| 18" | 1163011 | 1163211 | 1 |
| 24" | 1263011 | 1263211 | 1 |
| 30" | 1363011 | - | 2 |
| 36" | 1463011 | - | 2 |
| 42" | 1863011 | - | 2 |
| 48" | 1563011 | - | 2 |
| 513mm | 1663011 | - | 2 |
| 1000mm | 1763011 | - | 2 |

SUSPENDED COMPONENTS


Pencil Drawers and Pull Out Writing Tablets



| Width | Pencil Drawer | Pull Out Writing Tablet |
|-------|---------------|-------------------------|
| 15" | PEN0015 | PUL0015 |
| 18" | PEN0018 | PUL0018 |
| 20" | PEN0020 | PUL0020 |
| 24" | PEN0024 | PUL0024 |
| 30" | PEN0030 | PUL0030 |

- Pencil Drawers and Pull Out Writing Tablets suspends from tables with hanging rails. Both supplied with full extension runners.
- When ordering a table with aprons, add option code HU to suspend Pencil Drawer.

Hanging CPU Holder



| Item Number |
|-------------|
|-------------|

CPU1010

- Suspends from tables with hanging rails and is supplied with full extension runners to slide out CPU unit.
- Inside clear dimension of CPU holder is 9-3/16".
- Hanging CPU Holder does not fit on a 22" deep table.
- When ordering a table with aprons, add option code HU to suspend CPU holder.

SUSPENDED CABINET OPTIONS

| Option Code | Description |
|-------------|--|
| 00 | Modified door/drawer configuration |
| 01 | Full extension drawer runners |
| 02 | Label holders on all doors/drawers |
| 03 | Numeric or alpha identification disks on all doors/drawers (customer details required) |
| 04 | Standard louver configuration on doors (6 louvers at 3" wide x 5-1/2" long) |
| 05 | Non-standard louver configuration on doors (customer drawing required) |
| 07 | Cut-outs for electrical and/or plumbing services (customer drawing required) |
| 08 | Security proof panel |
| 09 | Stainless steel waste flap |
| 10 | One fixed drawer divider per drawer (centered in drawer) |
| 11 | Multiple drawer dividers per drawer (customer drawing required) |
| 12 | Pull out writing board - all cold rolled steel* |
| 13 | Pull out writing board - high pressure laminate writing surface* |
| 15 | Finished back on cabinet |
| 16 | Ship handles not installed |
| 18 | Horizontally mounted handles |
| 26 | 200 lbs full extension drawer runners |
| 40 | One extra cabinet shelf |
| 41 | Two extra cabinet shelves |
| 42 | One interior cabinet shelf on full extension runners (replaces fixed shelf) |
| 43 | Two interior cabinet shelves on full extension runners (replaces fixed shelf) |
| 44 | Three interior cabinet shelves on full extension runners (replaces fixed shelf) |
| 48 | Keyboard tray (replace drawer) |
| 50 | 304 stainless steel drawer bodies in powder coated cold rolled steel cabinet |
| 51 | Entire 304 stainless steel component |
| 52 | Entire 316 stainless steel component |
| 53 | Sigma Barrier™ galvanized steel coating |
| 54 | 304 stainless steel drawer bodies in Sigma Barrier™ galvanized steel cabinet |
| 56 | Other customer specified material |
| 82 | FRP liner for acid cabinets |
| 83 | Self closing and latching doors for solvent cabinet |
| 90 | Sink supports - for under counter mounted sinks |
| BL | "Base" label (use this option code to replace standard label on acid cabinets) |
| CL | Corrosive label |
| DC2 | Magnetic Door Catch |
| F1 | Welded finished top |
| F2 | Finished top and back |
| FCW | Front corners welded |
| FL | Fusible link with self closer for insulated solvent unit** |
| LA | Bilingual label (French/English) |
| LF | Longitudinal file |
| OL | Overlay design |
| PP | Pre-plumbed (option code not to be used for fume hoods)(customer details required) |
| PW | Pre-wired (option code not to be used for fume hoods) (customer details required) |
| PDC | Positive door catch |
| PFE | Pull out writing board on full extension runners* |
| RB | Removable back on drawer cabinet |
| SK | Skid item |
| SP | Silk screen printing (to replace standard labels on acid and/or flammable cabinets) |
| VK | Corrugated polyolefin vent kit for acid cabinets |
| VS | Vent holes on side of cabinet |
| T1 | Replace one shelf with one pull out tray on full extension runners with 2" high lip |
| T2 | Replace two shelves with two pull out trays on full extension runners with 2" high lip |
| T3 | Replace three shelves with three pull out trays on full extension runners with 2" high lip |
| WF | Wood fronts - see specifications |
| WO | Wood fronts provided and installed by others |

Notes: *30" wide cabinets and larger receive two pull out writing boards, one pull out board is provided on full width drawer cabinets.

Note: **Designed for use on solvent cabinets with self-closing and latching doors (option code 83). A fusible link holds doors open during use and melts at 165° F to automatically close doors in the event of a fire.