



SIGMA FLEX™

COUNTERTOP UPRIGHTS

Countertop Uprights

Countertop Uprights makes the most of the work surface space. Island or wall Countertop Uprights assemblies can be mounted to any new or existing work surface and suspend shelving components (accommodates a maximum of 6 Sigma Flex™ shelves). Uprights are 49" high and are offered in island or wall mounted assemblies in various widths.




Parts List

- 1) Top Shelf
- 2) Flex Post with welded mounting plate



Countertop Island Upright




Countertop Wall Upright

SIGMA FLEX™

49" HIGH COUNTERTOP UPRIGHTS


Island Countertop Uprights



Width	Item Number
24"	BIU1024
30"	BIU1030
36"	BIU1036
42"	BIU1042
48"	BIU1048
54"	BIU1054
60"	BIU1060
66"	BIU1066
72"	BIU1072

- 49" high. 12" deep.
- Countertop upright assemblies can be mounted to any new or existing work surface.
- Accommodates a maximum of 6 Sigma Flex™ shelves, adjustable in 1" increments.
- Ensure work surface is of sufficient strength to support additional load.
- Mounting requires drilling bolt holes in work surface (3/8" countersunk bolts recommended for mounting are included).

Island Countertop Upright Cover Plates



Island Upright Cover Panels	
With Service Cut-outs	Without Service Cut-outs
BIC0001	BIC0000

Island Upright Double Sided Cover Panels	
With Service Cut-outs	Without Service Cut-outs
BIC0003	BIC0002

- For BIU Island Countertop Uprights only.
- Used to cover visible screws that install a BIU.
- Available with and without service cut-out.

Wall Countertop Uprights



Width	Item Number
24"	BWU1024
30"	BWU1030
36"	BWU1036
42"	BWU1042
48"	BWU1048
54"	BWU1054
60"	BWU1060
66"	BWU1066
72"	BWU1072

- 49" high. 8-3/4" deep.
- Countertop upright assemblies can be mounted to any new or existing work surface.
- Accommodates a maximum of 3 Sigma Flex™ shelves, adjustable in 1" increments.
- Ensure work surface is of sufficient strength to support additional load.
- Mounting requires drilling bolt holes in work surface (3/8" countersunk bolts recommended for mounting are included).

Wall Countertop Upright Cover Plates



Item Number
BIC0004

- For BWU Wall Countertop Uprights only.
- Used to cover visible screws that install a BIU.