

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Open Units



| Width | Item Number |
|--------|-------------|
| 15" | 195000W |
| 18" | 115000W |
| 24" | 125000W |
| 30" | 135000W |
| 36" | 145000W |
| 42" | 185000W |
| 48" | 155000W |
| 513mm | 165000W |
| 1000mm | 175000W |

- Actual depths 7/8" less than nominal.

Full Door Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1950011W | 1950211W | 1 |
| 18" | 1150011W | 1150211W | 1 |
| 24" | 1250011W | 1250211W | 1 |
| 30" | 1350011W | - | 2 |
| 36" | 1450011W | - | 2 |
| 42" | 1850011W | - | 2 |
| 48" | 1550011W | - | 2 |
| 513mm | 1650011W | 1650211W | 1 |
| 1000mm | 1750011W | - | 2 |

- Cabinets receive one shelf.

Door/Drawer Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1950022W | 1950222W | 1 |
| 18" | 1150022W | 1150222W | 1 |
| 24" | 1250022W | 1250222W | 1 |
| 30" | 1350022W | - | 2 |
| 36" | 1450022W | - | 2 |
| 42" | 1850022W | - | 2 |
| 48" | 1550022W | - | 2 |
| 513mm | 1650022W | 1650222W | 1 |
| 1000mm | 1750022W | - | 2 |

- Cabinets up to 24" wide receive one shelf.
- Cabinets 25" and wider do not have a shelf.



| Width | Right Hand | Left Hand |
|--------|------------|-----------|
| 30" | 1350041W | 1350014W |
| 36" | 1450041W | 1450014W |
| 42" | 1850041W | 1850014W |
| 48" | 1550041W | 1550014W |
| 1000mm | 1750041W | 1750014 W |

- Left hand Door/Drawer Unit shown.
- Right hand Door/Drawer Units are reverse of what is shown.
- Cabinets receive one shelf.

Sliding Door Units



| Width | Item Number |
|--------|-------------|
| 30" | 1350711W |
| 36" | 1450711W |
| 42" | 1850711W |
| 48" | 1550711W |
| 1000mm | 1750711W |

- Actual depths 7/8" less than nominal.
- Cabinets receive one shelf.

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950044W |
| 18" | 1150044W |
| 24" | 1250044W |
| 30" | 1350044W |
| 36" | 1450044W |
| 42" | 1850044W |
| 48" | 1550044W |
| 513mm | 1650044W |
| 1000mm | 1750044W |



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950055W |
| 18" | 1150055W |
| 24" | 1250055W |
| 30" | 1350055W |
| 36" | 1450055W |
| 42" | 1850055W |
| 48" | 1550055W |
| 513mm | 1650055W |
| 1000mm | 1750055W |



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950066W |
| 18" | 1150066W |
| 24" | 1250066W |
| 30" | 1350066W |
| 36" | 1450066W |
| 42" | 1850066W |
| 48" | 1550066W |
| 513mm | 1650066W |
| 1000mm | 1750066W |



| Width | Item Number |
|--------|-------------|
| 30" | 1350097W |
| 36" | 1450097W |
| 42" | 1850097W |
| 48" | 1550097W |
| 1000mm | 1750097W |



| Width | Item Number |
|--------|-------------|
| 30" | 1350098W |
| 36" | 1450098W |
| 42" | 1850098W |
| 48" | 1550098W |
| 1000mm | 1750098W |

SUSPENDED CABINETS

24" HIGH X 22" DEEP

File Drawer Units



up to 24" wide



30" and wider



up to 24" wide



30" and wider



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1950122W |
| 18" | 1150122W |
| 24" | 1250122W |
| 30" | 1350122W |
| 36" | 1450122W |
| 42" | 1850122W |
| 48" | 1550122W |
| 513mm | 1650122W |
| 1000mm | 1750122W |

| Width | Item Number |
|--------|-------------|
| 15" | 1950177W |
| 18" | 1150177W |
| 24" | 1250177W |
| 30" | 1350177W |
| 36" | 1450177W |
| 42" | 1850177W |
| 48" | 1550177W |
| 513mm | 1650177W |
| 1000mm | 1750177W |

| Width | Item Number |
|--------|-------------|
| 15" | 1950188W |
| 18" | 1150188W |
| 24" | 1250188W |
| 30" | 1350188W |
| 36" | 1450188W |
| 42" | 1850188W |
| 48" | 1550188W |
| 513mm | 1650188W |
| 1000mm | 1750188W |

Sink Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1953011W | 1953211W | 1 |
| 18" | 1153011W | 1153211W | 1 |
| 24" | 1253011W | 1253211W | 1 |
| 30" | 1353011W | - | 2 |
| 36" | 1453011W | - | 2 |
| 42" | 1853011W | - | 2 |
| 48" | 1553011W | - | 2 |
| 513mm | 1653011W | 1653211W | 1 |
| 1000mm | 1753011W | - | 2 |

• Add option code 90 to order Sink Hanger Supports.

Acid Storage Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1951011W | 1951211W | 1 |
| 18" | 1151011W | 1151211W | 1 |
| 24" | 1251011W | 1251211W | 1 |
| 30" | 1351011W | - | 2 |
| 36" | 1451011W | - | 2 |
| 42" | 1851011W | - | 2 |
| 48" | 1551011W | - | 2 |
| 513mm | 1651011W | 1651211W | 1 |
| 1000mm | 1751011W | - | 2 |

- Cabinet receives a molded liner and one chemically resistant shelf.
- Add option code VK to order a Corrugated Polyolefin Vent Kit for these units.
- Note: Acid Cabinets must be vented. Mott recommends venting acid cabinets with a Polyolefin Vent Kit (PVSA for FRP lined cabinets).

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Solvent Storage Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1952011W | 1952211W | 1 |
| 18" | 1152011W | 1152211W | 1 |
| 24" | 1252011W | 1252211W | 1 |
| 30" | 1352011W | - | 2 |
| 36" | 1452011W | - | 2 |
| 42" | 1852011W | - | 2 |
| 48" | 1552011W | - | 2 |
| 513mm | 1652011W | 1652211W | 1 |
| 1000mm | 1752011W | - | 2 |

- Solvent Units are designed to meet OSHA and NFPA 30 standards.
- Add option code 83 to order Self Closing/Latching doors.
- Add option code "FL" to order Self Closing/latching doors with a fusible link.
- See page L10 for Solvent cabinet details.

SUSPENDED CABINETS

18" HIGH X 22" DEEP

Open Units



| Width | Item Number |
|--------|-------------|
| 15" | 1960000W |
| 18" | 1160000W |
| 24" | 1260000W |
| 30" | 1360000W |
| 36" | 1460000W |
| 42" | 1860000W |
| 48" | 1560000W |
| 513mm | 1660000W |
| 1000mm | 1760000W |

- Actual depths 7/8" less than nominal.

Full Door Unit



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1960011W | 1960211W | 1 |
| 18" | 1160011W | 1160211W | 1 |
| 24" | 1260011W | 1260211W | 1 |
| 30" | 1360011W | - | 2 |
| 36" | 1460011W | - | 2 |
| 42" | 1860011W | - | 2 |
| 48" | 1560011W | - | 2 |
| 513mm | 1660011W | 1660211W | 1 |
| 1000mm | 1760011W | - | 2 |

Door/Drawer Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1960022W | 1960222W | 1 |
| 18" | 1160022W | 1160222W | 1 |
| 24" | 1260022W | 1260222W | 1 |
| 30" | 1360022W | - | 2 |
| 36" | 1460022W | - | 2 |
| 42" | 1860022W | - | 2 |
| 48" | 1560022W | - | 2 |
| 513mm | 1660022W | 1660222W | 1 |
| 1000mm | 1760022W | - | 2 |



| Width | Right Hinged | Left Hinged |
|--------|--------------|-------------|
| 30" | 1360041W | 1360014W |
| 36" | 1460041W | 1460014W |
| 42" | 1860041W | 1860014W |
| 48" | 1560041W | 1560014W |
| 1000mm | 1760041W | 1760014W |

- Right hand Door/Drawer Unit shown. Left hand reversed.

SUSPENDED CABINETS

18" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1960044W |
| 18" | 1160044W |
| 24" | 1260044W |
| 30" | 1360044W |
| 36" | 1460044W |
| 42" | 1860044W |
| 48" | 1560044W |
| 513mm | 1660044W |
| 1000mm | 1760044W |

File Drawer Units



up to 24" wide



30" and wider

| Width | Item Number |
|--------|-------------|
| 15" | 1960177W |
| 18" | 1160177W |
| 24" | 1260177W |
| 30" | 1360177W |
| 36" | 1460177W |
| 42" | 1860177W |
| 48" | 1560177W |
| 513mm | 1660177W |
| 1000mm | 1760177W |

Sink Units



up to 24" wide



30" and wider

| Width | Right Hinged | Left Hinged | # of Doors |
|--------|--------------|-------------|------------|
| 15" | 1963011W | 1963211W | 1 |
| 18" | 1163011W | 1163211W | 1 |
| 24" | 1263011W | 1263211W | 1 |
| 30" | 1363011W | - | 2 |
| 36" | 1463011W | - | 2 |
| 42" | 1863011W | - | 2 |
| 48" | 1563011W | - | 2 |
| 513mm | 1663011W | 1663211W | 1 |
| 1000mm | 1763011W | - | 2 |

- Add option code 90 to order Sink Hanger Supports.

Hanging Rails



| Width | Item Number |
|--------|-------------|
| 15" | HR00015 |
| 18" | HR00018 |
| 24" | HR00024 |
| 30" | HR00030 |
| 36" | HR00036 |
| 42" | HR00042 |
| 48" | HR00048 |
| 513mm | HR00020 |
| 1000mm | HR00040 |

- Pair of hanging rails only suitable for wood suspended cabinets.

SUSPENDED CABINET OPTIONS

| Option Code | Description |
|-------------|--|
| 00 | Modified door/drawer configuration |
| 02 | Label holders on all doors/drawers |
| 03 | Numeric or alpha identification disks on all doors/drawers (customer details required) |
| 04 | Standard louver configuration on doors (6 louvers at 3" wide x 5-1/2" long) |
| 05 | Non-standard louver configuration on doors (customer drawing required) |
| 07 | Cut-outs for electrical and/or plumbing services (customer drawing required) |
| 08 | Security proof panel |
| 09 | Stainless steel waste flap |
| 10 | One fixed drawer divider per drawer (centered in drawer) |
| 11 | Multiple drawer dividers per drawer (customer drawing required) |
| 12 | Pull out writing board - all cold rolled steel* |
| 13 | Pull out writing board - high pressure laminate writing surface* |
| 15 | Finished back on cabinet |
| 16 | Ship handles not installed |
| 18 | Horizontally mounted handles |
| 26 | 200 lbs full extension drawer runners |
| 40 | One extra cabinet shelf |
| 41 | Two extra cabinet shelves |
| 42 | One interior cabinet shelf on full extension runners (replaces fixed shelf) |
| 43 | Two interior cabinet shelves on full extension runners (replaces fixed shelf) |
| 44 | Three interior cabinet shelves on full extension runners (replaces fixed shelf) |
| 48 | Keyboard tray (replace drawer) |
| 82 | FRP liner for acid cabinets |
| 83 | Self closing and latching doors for solvent cabinet |
| 90 | Sink supports - for under counter mounted sinks |
| BL | "Base" label (use this option code to replace standard label on acid cabinets) |
| CL | Corrosive label |
| DC2 | Magnetic Door Catch |
| F2 | Finished top and back |
| FL | Fusible link with self closer for insulated unit** |
| LA | Bilingual label (French/English) |
| LF | Longitudinal file |
| OLD | Actual depth |
| PP | Pre-plumbed (option code not to be used for fume hoods)(customer details required) |
| PW | Pre-wired (option code not to be used for fume hoods) (customer details required) |
| PDC | Positive door catch |
| PFE | Pull out writing board on full extension runners* |
| RB | Removable back on drawer cabinet |
| SK | Skid item |
| SP | Silk screen printing (to replace standard labels on acid and/or flammable cabinets) |
| VK | Corrugated polyolefin vent kit for acid cabinets |
| VS | Vent holes on side of cabinet |
| T1 | Replace one shelf with one pull out tray on full extension runners with 2" high lip |
| T2 | Replace two shelves with two pull out trays on full extension runners with 2" high lip |
| T3 | Replace three shelves with three pull out trays on full extension runners with 2" high lip |

Notes: *30" wide cabinets and above receive 2 pull out writing boards unless it's a full width drawer cabinet.

Note: **Designed for use on solvent cabinets with self-closing and latching doors (option code 83). A fusible link holds doors open during use and melts at 165° F to automatically close doors in the event of a fire.