

MOBILE CABINETS & TABLES

MOBILE FURNITURE SOLUTIONS

The Solution to Change...


Laboratory technology changes today at an ever-increasing rate. Today's challenge is to provide the longest possible life to your lab. Mott's flexible furniture systems successfully meet this challenge by delivering ultimate adaptability and adjustability with exceptional architectural style.

Mobile Cabinets



Mobile cabinets are available in numerous configurations and various sizes. Cabinets are available in powder coated steel, wood, wood fronts on powder coated steel or stainless steel in either contemporary inset or architectural overlay design. Heavy duty locking 3-1/2" twin wheel casters provide a total load rating (including mobile cabinet and work top) of 500lbs. Tip resisting counter weights are provided to aid in safe operation. Resulting mobility allows the user to quickly and easily move to another workstation, as tasks change.

Mobile Cabinet Height With Casters



Mobile Tables



Mobile tables are ideal work surface solutions where overhead storage or service carriers already exist, or when work surface mobility is a critical requirement. Offered in numerous configurations and sizes, tables are available in powder coated steel, stainless steel or wood and can be configured to include lower suspended cabinets and shelving. Ideal for sitting or standing activities, these can be adjustable to be readily reconfigured for changing applications or process technologies.

For additional information and help in planning contact Mott Manufacturing.