

**MISCELLANEOUS ITEMS**

**Sigma Frame™ Shelf Retainer Rods**

	Steel Shelves			Wood Shelves		Retainer Rod Parts
	Width	Powder Coated	Stainless Steel	Powder Coated	Stainless Steel	
 ≤ 30" wide assemblies	24"	PRR0024-MI	SRRO024-MI	PRR0024-MIW	SRRO024-MIW	1 rod & 2 turrets
	30"	PRR0030-MI	SRRO030-MI	PRR0030-MIW	SRRO030-MIW	1 rod & 2 turrets
	36"	PRR0036-MI	SRRO036-MI	PRR0036-MIW	SRRO036-MIW	1 rod & 3 turrets
 > 30" wide assemblies	42"	PRR0042-MI	SRRO042-MI	PRR0042-MIW	SRRO042-MIW	1 rod & 3 turrets
	48"	PRR0048-MI	SRRO048-MI	PRR0048-MIW	SRRO048-MIW	1 rod & 3 turrets
	60"	PRR0060-MI	SRRO060-MI	PRR0060-MIW	SRRO060-MIW	1 rod & 3 turrets
	513mm	PRR000A-MI	SRRO00A-MI	PRR000A-MIW	SRRO00A-MIW	1 rod & 2 turrets
	1000mm	PRR000B-MI	SRRO00B-MI	PRR000B-MIW	SRRO00B-MIW	1 rod & 3 turrets
	1512mm	PRR000C-MI	SRRO00C-MI	PRR000C-MIW	SRRO00C-MIW	1 rod & 3 turrets
2000mm	PRR000D-MI	SRRO00D-MI	PRR000D-MIW	SRRO00D-MIW	1 rod & 4 turrets	

- 1" high retainer rod (5/16" dia.) Prevents articles from falling off the shelf.
- Available in stainless steel or powder coated stainless steel and can be added to existing shelves.
- Threaded turrets includes hardware for mounting to shelf (shelf must have holes drilled to accept turrets).
- Ships loose.
- MI and MIW at the end of the part numbers represent option codes.


**Sigma Frame™ Shelf Brackets For Wall Pilasters**



Depth	Standard Profile	Radius Profile	Tapered Profile
10"	SBP0010	SBP3010	SBP4010
12"	SBP0012	SBP3012	SBP4012
14"	SBP0014	SBP3014	SBP4014
16"	SBP0016	SBP3016	SBP4016
18"	SBP0018	SBP3018	SBP4018

- Brackets only suitable for Sigma Frame™ and "WPA" components.
- Pair of shelf brackets are 6" high.
- Add option code NN for no alignment tabs on SBP brackets.

**Sigma Frame™ Inverted Shelf Brackets For Wall Pilasters**



- Brackets only suitable for Reagent components.
- Pair of shelf brackets are 4-1/4" high.

Depth	Tapered Profile	Radius Profile
10"	SBP1010	SBP1110
12"	SBP1012	SBP1112
14"	SBP1014	SBP1114
16"	SBP1016	SBP1116
18"	SBP1018	SBP1118

**Under Mount Shelf Task Lights**



- Order light a minimum 6" shorter than the shelf.
- Low profile (.98" wide by 1.65" high) T5 fluorescent task light mounts to underside of the shelf.
- Shatterproof polycarbonate lens protects the lamp.
- Easily linkable in one continuous line using Flexible Connectors. 6" Flexible Connector comes standard with each light. See below to order longer Flexible Connectors.
- On-off rocker switch located on the side of the fixture.
- 8000 hour, cool white replaceable T5 lamp included.
- UL listed.
- Ships loose, includes Mott's mounting brackets.

Length of Light	Item Number
14"	TLW0014
23"	TLW0023
36"	TLW0036
48"	TLW0048
59"	TLW0059

**Task Light Flexible Connectors**



- Flexible Connector allows to easily link T5 fluorescent lights in one continuous line.
- Ships loose.

Length	Item Number
12"	LFC0012
24"	LFC0024
36"	LFC0036