

BASE CABINETS

34-3/4" HIGH X 22" DEEP

Open Units



Width	Item Number
12"	1010000
15"	1910000
18"	1110000
24"	1210000
30"	1310000
36"	1410000
42"	1810000
48"	1510000
513mm	1610000
1000mm	1710000

Full Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1010011	1010211	1
15"	1910011	1910211	1
18"	1110011	1110211	1
24"	1210011	1210211	1
30"	1310011	-	2
36"	1410011	-	2
42"	1810011	-	2
48"	1510011	-	2
513mm	1610011	1610211	1
1000mm	1710011	-	2

Door/Drawer Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1010022	1010222	1
15"	1910022	1910222	1
18"	1110022	1110222	1
24"	1210022	1210222	1
30"	1310022	-	2
36"	1410022	-	2
42"	1810022	-	2
48"	1510022	-	2
513mm	1610022	1610222	1
1000mm	1710022	-	2



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1010033	1010233	1
15"	1910033	1910233	1
18"	1110033	1110233	1
24"	1210033	1210233	1
30"	1310033	-	2
36"	1410033	-	2
42"	1810033	-	2
48"	1510033	-	2
513mm	1610033	1610233	1
1000mm	1710033	-	2



Width	Item Number
30"	1310093
36"	1410093
42"	1810093
48"	1510093
1000mm	1710093



Width	Item Number
30"	1310096
36"	1410096
42"	1810096
48"	1510096
1000mm	1710096

BASE CABINETS

34-3/4" HIGH X 22" DEEP

Door/Drawer Units



Width	Right Hand	Left Hand
30"	1310032	1310023
36"	1410032	1410023
42"	1810032	1810023
48"	1510032	1510023
1000mm	1710032	1710023

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310021	1310012
36"	1410021	1410012
42"	1810021	1810012
48"	1510021	1510012
1000mm	1710021	1710012

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310031	1310013
36"	1410031	1410013
42"	1810031	1810013
48"	1510031	1510013
1000mm	1710031	1710013

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310041	1310014
36"	1410041	1410014
42"	1810041	1810014
48"	1510041	1510014
1000mm	1710041	1710014

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310051	1310015
36"	1410051	1410015
42"	1810051	1810015
48"	1510051	1510015
1000mm	1710051	1710015

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310061	1310016
36"	1410061	1410016
42"	1810061	1810016
48"	1510061	1510016
1000mm	1710061	1710016

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310071	1310017
36"	1410071	1410017
42"	1810071	1810017
48"	1510071	1510017
1000mm	1710071	1710017

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310081	1310018
36"	1410081	1410018
42"	1810081	1810018
48"	1510081	1510018
1000mm	1710081	1710018

- Left hand door/drawer unit shown
- Right hand reversed
- If locks are required, the lowest drawer of the cabinet will require a slot in the bottom of the drawer body to accommodate the lock cam



Width	Right Hand	Left Hand
30"	13100B1	131001B
36"	14100B1	141001B
42"	18100B1	181001B
48"	15100B1	151001B
1000mm	17100B1	171001B

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310942	1310924
36"	1410942	1410924
42"	1810942	1810924
48"	1510942	1510924
1000mm	1710942	1710924

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310042	1310024
36"	1410042	1410024
42"	1810042	1810024
48"	1510042	1510024
1000mm	1710042	1710024

- Left hand door/drawer unit shown
- Right hand reversed

BASE CABINETS
34-3/4" HIGH X 22" DEEP
Door/Drawer Units


Width	Right Hand	Left Hand
30"	1310052	1310025
36"	1410052	1410025
42"	1810052	1810025
48"	1510052	1510025
1000mm	1710052	1710025

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310062	1310026
36"	1410062	1410026
42"	1810062	1810026
48"	1510062	1510026
1000mm	1710062	1710026

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310072	1310027
36"	1410072	1410027
42"	1810072	1810027
48"	1510072	1510027
1000mm	1710072	1710027

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310082	1310028
36"	1410082	1410028
42"	1810082	1810028
48"	1510082	1510028
1000mm	1710082	1710028

- Left hand door/drawer unit shown
- Right hand reversed
- If locks are required, the lowest drawer of the cabinet will require a slot in the bottom of the drawer body to accommodate the lock cam



Width	Right Hand	Left Hand
30"	1310043	1310034
36"	1410043	1410034
42"	1810043	1810034
48"	1510043	1510034
1000mm	1710043	1710034

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310063	1310036
36"	1410063	1410036
42"	1810063	1810036
48"	1510063	1510036
1000mm	1710063	1710036

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310053	1310035
36"	1410053	1410035
42"	1810053	1810035
48"	1510053	1510035
1000mm	1710053	1710035

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1310083	1310038
36"	1410083	1410038
42"	1810083	1810038
48"	1510083	1510038
1000mm	1710083	1710038

- Left hand door/drawer unit shown
- Right hand reversed
- If locks are required, the lowest drawer of the cabinet will require a slot in the bottom of the drawer body to accommodate the lock cam



Width	Right Hand	Left Hand
30"	1310073	1310037
36"	1410073	1410037
42"	1810073	1810037
48"	1510073	1510037
1000mm	1710073	1710037

- Left hand door/drawer unit shown
- Right hand reversed

BASE CABINETS

34-3/4" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider



up to 24" wide



30" and wider



up to 24" wide



30" and wider



up to 24" wide



30" and wider



up to 24" wide



30" and wider

Width	Item Number
12"	1010090
15"	1910090
18"	1110090
24"	1210090
30"	13100NN
36"	14100NN
42"	18100NN
48"	15100NN
513mm	1610090
1000Mm	17100NN

Width	Item Number
12"	1010044
15"	1910044
18"	1110044
24"	1210044
30"	1310044
36"	1410044
42"	1810044
48"	1510044
513mm	1610044
1000mm	1710044

Width	Item Number
15"	19100BB
18"	11100BB
24"	12100BB
30"	13100BB
36"	14100BB
42"	18100BB
48"	15100BB
513mm	16100BB
1000mm	17100BB

Width	Item Number
15"	1910055
18"	1110055
24"	1210055
30"	1310055
36"	1410055
42"	1810055
48"	1510055
513mm	1610055
1000mm	1710055

Width	Item Number
15"	1910066
18"	1110066
24"	1210066
30"	1310066
36"	1410066
42"	1810066
48"	1510066
513mm	1610066
1000mm	1710066

BASE CABINETS

34-3/4" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider

Width	Item Number
15"	1910077
18"	1110077
24"	1210077
30"	1310077
36"	1410077
42"	1810077
48"	1510077
513mm	1610077
1000mm	1710077



up to 24" wide



30" and wider

Width	Item Number
15"	1910088
18"	1110088
24"	1210088
30"	1310088
36"	1410088
42"	1810088
48"	1510088
513mm	1610088
1000mm	1710088

- If locks are required, the lowest drawer of the cabinet will require a slot in the bottom of the drawer body to accommodate the lock cam



Width	Item Number
30"	1310A91
36"	1410A91
42"	1810A91
48"	1510A91
1000mm	1710A91



Width	Item Number
30"	1310099
36"	1410099
42"	1810099
48"	1510099
1000mm	1710099



Width	Item Number
30"	1310090
36"	1410090
42"	1810090
48"	1510090
1000mm	1710090



Width	Item Number
30"	1310B44
36"	1410B44
42"	1810B44
48"	1510B44
1000mm	1710B44



Width	Item Number
30"	1310097
36"	1410097
42"	1810097
48"	1510097
1000mm	1710097



Width	Right Hand	Left Hand
30"	1310094	1310092
36"	1410094	1410092
42"	1810094	1810092
48"	1510094	1510092
1000mm	1710094	1710092

- Left hand drawer unit shown
- Right hand reversed

BASE CABINETS

34-3/4" HIGH X 22" DEEP

Drawer Units



Width	Item Number
30"	1310098
36"	1410098
42"	1810098
48"	1510098
1000mm	1710098



Width	Item Number
30"	1310091
36"	1410091
42"	1810091
48"	1510091
1000mm	1710091



Width	Item Number
30"	1310039
36"	1410039
42"	1810039
48"	1510039
1000mm	1710039



Width	Item Number
30"	1310095
36"	1410095
42"	1810095
48"	1510095
1000mm	1710095



Width	Item Number
30"	13100CC
36"	14100CC
42"	18100CC
48"	15100CC
1000mm	17100CC



Width	Right Hand	Left Hand
30"	1310019	1310029
36"	1410019	1410029
42"	1810019	1810029
48"	1510019	1510029
1000mm	1710019	1710029

- Left hand drawer unit shown
- Right hand reversed

File Drawer Units



up to 24" wide



30" and wider



up to 24" wide



30" and wider

Width	Item Number
12"	1010133
15"	1910133
18"	1110133
24"	1210133
30"	1310133
36"	1410133
42"	1810133
48"	1510133
513mm	1610133
1000mm	1710133

Width	Item Number
12"	1010100
15"	1910100
18"	1110100
24"	1210100
30"	1310100
36"	1410100
42"	1810100
48"	1510100
513mm	1610100
1000mm	1710100

BASE CABINETS

34-3/4" HIGH X 22" DEEP

File Drawer Units



up to 24" wide



30" and wider

Width	Item Number
12"	1010111
15"	1910111
18"	1110111
24"	1210111
30"	1310111
36"	1410111
42"	1810111
48"	1510111
513mm	1610111
1000mm	1710111



up to 24" wide



30" and wider

Width	Item Number
12"	1010122
15"	1910122
18"	1110122
24"	1210122
30"	1310122
36"	1410122
42"	1810122
48"	1510122
513mm	1610122
1000mm	1710122



up to 24" wide



30" and wider

Width	Item Number
12"	10101AA
15"	19101AA
18"	11101AA
24"	12101AA
30"	13101AA
36"	14101AA
42"	18101AA
48"	15101AA
513mm	16101AA
1000mm	17101AA

BASE CABINETS

34-3/4" HIGH X 22" DEEP

Waste/Storage Bin Units



up to 24" wide



30" and wider



Tilt Out Unit
(shown with
optional bin)




Pull Out Unit

Width	Tilt Out Units	Pull Out Units	Bin Capacity
15"	1910500	1910511	1
18"	1110500	1110511	1
24"	1210500	1210511	1
30"	1310500	1310511	2
36"	1410500	1410511	2
42"	1810500	1810511	2
48"	1510500	1510511	2
513mm	1610500	1610511	1
1000mm	1710500	1710511	2

- Add option code 84 to order stainless steel waste/storage bin, or add option code 85 for PVC waste/storage bin for Tilt Out units
- Tilt out units are hinged at the bottom
- Pull out units include full extension runners

Diagonal Corner Units



Description	Right Hinged	Left Hinged
Full Door Unit	1017011	1017111
False Panel Unit	1017422	1017922


- Left hinged, full door corner unit shown
- Add option code 89 to order a Lazy Susan

Sliding Door Units



Width	Item Number
30"	1310711
36"	1410711
42"	1810711
48"	1510711
1000mm	1710711

Blind Corner Units



Width	Door on Right	Door on Left
30"	1310311	1310411
36"	1410311	1410411
42"	1810311	1810411
48"	1510311	1510411
1000mm	1710311	1710411

- Right door blind corner unit shown, left door reversed
- Center stile located at the middle point

BASE CABINETS
34-3/4" HIGH X 22" DEEP
Sink Units


up to 24" wide



30" and wider



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1013011	1013211	1
15"	1913011	1913211	1
18"	1113011	1113211	1
24"	1213011	1213211	1
30"	1313011	-	2
36"	1413011	-	2
42"	1813011	-	2
48"	1513011	-	2
513mm	1613011	1613211	1
1000mm	1713011	-	2


Width	Right Hinged	Left Hinged	# of Doors
12"	1018022	1018222	1
15"	1918022	1918222	1
18"	1118022	1118222	1
24"	1218022	1218222	1
30"	1318022	-	2
36"	1418022	-	2
42"	1818022	-	2
48"	1518022	-	2
513mm	1618022	1618222	1
1000mm	1718022	-	2

- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

End Sink Units


Description	54" Wide	58" Wide
Full Door Unit	9900017	9900018
False Panel Unit	9900021	9900022

- Full door end sink unit shown
- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

Diagonal Corner Sink Units


Top View

Description	Right Hinged	Left Hinged
Full Door Unit	1017211	1017511
False Panel Unit	1017722	1017622

- Left hinged, full door corner sink unit shown
- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

BASE CABINETS

28-3/4" HIGH X 22" DEEP

Open Units



Width	Item Number
12"	1020000
15"	1920000
18"	1120000
24"	1220000
30"	1320000
36"	1420000
42"	1820000
48"	1520000
513mm	1620000
1000mm	1720000

Full Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1020011	1020211	1
15"	1920011	1920211	1
18"	1120011	1120211	1
24"	1220011	1220211	1
30"	1320011	-	2
36"	1420011	-	2
42"	1820011	-	2
48"	1520011	-	2
513mm	1620011	1620211	1
1000mm	1720011	-	2

Door/Drawer Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1020022	1020222	1
15"	1920022	1920222	1
18"	1120022	1120222	1
24"	1220022	1220222	1
30"	1320022	-	2
36"	1420022	-	2
42"	1820022	-	2
48"	1520022	-	2
513mm	1620022	1620222	1
1000mm	1720022	-	2



up to 24" wide



30" and wider

Width	Item Number	Left Hinged	# of Doors
12"	1020033	1020233	1
15"	1920033	1920233	1
18"	1120033	1120233	1
24"	1220033	1220233	1
30"	1320033	-	2
36"	1420033	-	2
42"	1820033	-	2
48"	1520033	-	2
513mm	1620033	1620233	1
1000mm	1720033	-	2



Width	Item Number
30"	1320093
36"	1420093
42"	1820093
48"	1520093
1000mm	1720093



Width	Item Number
30"	1320096
36"	1420096
42"	1820096
48"	1520096
1000mm	1720096

BASE CABINETS
28-3/4" HIGH X 22" DEEP
Door/Drawer Units


Width	Right Hand	Left Hand
30"	1320032	1320023
36"	1420032	1420023
42"	1820032	1820023
48"	1520032	1520023
1000mm	1720032	1720023

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320942	1320924
36"	1420942	1420924
42"	1820942	1820924
48"	1520942	1520924
1000mm	1780942	1720924

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320021	1320012
36"	1420021	1420012
42"	1820021	1820012
48"	1520021	1520012
1000mm	1720021	1720012

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320042	1320024
36"	1420042	1420024
42"	1820042	1820024
48"	1520042	1520024
1000mm	1720042	1720024

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320031	1320013
36"	1420031	1420013
42"	1820031	1820013
48"	1520031	1520013
1000mm	1720031	1720013

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320052	1320025
36"	1420052	1420025
42"	1820052	1820025
48"	1520052	1520025
1000mm	1720052	1720025

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320041	1320014
36"	1420041	1420014
42"	1820041	1820014
48"	1520041	1520014
1000mm	1720041	1720014

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320062	1320026
36"	1420062	1420026
42"	1820062	1820026
48"	1520062	1520026
1000mm	1720062	1720026

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320051	1320015
36"	1420051	1420015
42"	1820051	1820015
48"	1520051	1520015
1000mm	1720051	1720015

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320072	1320027
36"	1420072	1420027
42"	1820072	1820027
48"	1520072	1520027
1000mm	1720072	1720027

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320061	1320016
36"	1420061	1420016
42"	1820061	1820016
48"	1520061	1520016
1000mm	1720061	1720016

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320043	1320034
36"	1420043	1420034
42"	1820043	1820034
48"	1520043	1520034
1000mm	1720043	1720034

- Left hand door/drawer unit shown
- Right hand reversed

BASE CABINETS

28-3/4" HIGH X 22" DEEP

Door/Drawer Units



Width	Right Hand	Left Hand
30"	1320053	1320035
36"	1420053	1420035
42"	1820053	1820035
48"	1520053	1520035
1000mm	1720053	1720035

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320063	1320036
36"	1420063	1420036
42"	1820063	1820036
48"	1520063	1520036
1000mm	1720063	1720036

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1320073	1320037
36"	1420073	1420037
42"	1820073	1820037
48"	1520073	1520037
1000mm	1720073	1720037

- Left hand door/drawer unit shown
- Right hand reversed

Drawer Units



up to 24" wide



30" and wider



up to 24" wide



30" and wider



up to 24" wide



30" and wider

Width	Item Number
12"	1020044
15"	1920044
18"	1120044
24"	1220044
30"	1320044
36"	1420044
42"	1820044
48"	1520044
513mm	1620044
1000mm	1720044

Width	Item Number
12"	10200MM
15"	19200MM
18"	11200MM
24"	12200MM
30"	13200MM
36"	14200MM
42"	18200MM
48"	15200MM
513mm	16200MM
1000mm	17200MM

Width	Item Number
12"	1020055
15"	1920055
18"	1120055
24"	1220055
30"	1320055
36"	1420055
42"	1820055
48"	1520055
513mm	1620055
1000mm	1720055

BASE CABINETS

28-3/4" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider

Width	Item Number
15"	1920066
18"	1120066
24"	1220066
30"	1320066
36"	1420066
42"	1820066
48"	1520066
513mm	1620066
1000mm	1720066



up to 24" wide



30" and wider

Width	Item Number
15"	1920077
18"	1120077
24"	1220077
30"	1320077
36"	1420077
42"	1820077
48"	1520077
513mm	1620077
1000mm	1720077



Width	Item Number
30"	1320097
36"	1420097
42"	1820097
48"	1520097
1000mm	1720097



Width	Item Number
30"	13200LL
36"	14200LL
42"	18200LL
48"	15200LL
1000mm	17200LL



Width	Item Number
30"	1320098
36"	1420098
42"	1820098
48"	1520098
1000mm	1720098



Width	Item Number
30"	1320091
36"	1420091
42"	1820091
48"	1520091
1000mm	1720091



Width	Item Number
30"	1320099
36"	1420099
42"	1820099
48"	1520099
1000mm	1720099



Width	Right Hand	Left Hand
30"	1320094	1320092
36"	1420094	1420092
42"	1820094	1820092
48"	1520094	1520092
1000mm	1720094	1720092

- Left hand door/drawer unit shown
- Right hand reversed

BASE CABINETS

28-3/4" HIGH X 22" DEEP

Drawer Units



Width	Item Number
30"	1320095
36"	1420095
42"	1820095
48"	1520095
1000mm	1720095



Width	Item Number
30"	1320B44
36"	1420B44
42"	1820B44
48"	1520B44
1000mm	1720B44

File Drawer Units



up to 24" wide



30" and wider

Width	Item Number
12"	1020122
15"	1920122
18"	1120122
24"	1220122
30"	1320122
36"	1420122
42"	1820122
48"	1520122
513mm	1620122
1000mm	1720122



up to 24" wide



30" and wider

Width	Item Number
12"	1020177
15"	1920177
18"	1120177
24"	1220177
30"	1320177
36"	1420177
42"	1820177
48"	1520177
513mm	1620177
1000mm	1720177



up to 24" wide



30" and wider

Width	Item Number
12"	1020188
15"	1920188
18"	1120188
24"	1220188
30"	1320188
36"	1420188
42"	1820188
48"	1520188
513mm	1620188
1000mm	1720188

BASE CABINETS

28-3/4" HIGH X 22" DEEP

Waste/Storage Bin Units



up to 24" wide



30" and wider



Tilt Out Unit
(shown with
optional bin)



Pull Out Unit

Width	Tilt Out Units	Pull Out Units	Bin Capacity
15"	1920500	1920511	1
18"	1120500	1120511	1
24"	1220500	1220511	1
30"	1320500	1320511	2
36"	1420500	1420511	2
42"	1820500	1820511	2
48"	1520500	1520511	2
513mm	1620500	1620511	1
1000mm	1720500	1720511	2


- Add option code 84 to order stainless steel waste/storage bin, or add option code 85 for PVC waste/storage bin for Tilt Out units
- Tilt out units are hinged at the bottom
- Pull out units include full extension runners

Sliding Door Units



Width	Item Number
30"	1320711
36"	1420711
42"	1820711
48"	1520711
1000mm	1720711


Blind Corner Units



Width	Door on Right	Door on Left
30"	1320311	1320411
36"	1420311	1420411
42"	1820311	1820411
48"	1520311	1520411
1000mm	1720311	1720411

- Right door blind corner unit shown, left door reversed
- Center stile located at the middle point

Diagonal Corner Units



Description	Right Hinged	Left Hinged
Full Door Unit	1027011	1027111
False Panel Unit	1027422	1027922

- Left hinged, full door corner unit shown
- Add option code 89 to order a Lazy Susan

BASE CABINETS

28-3/4" HIGH X 22" DEEP

Sink Units



up to 24" wide



30" and wider



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1023011	1023211	1
15"	1923011	1923211	1
18"	1123011	1123211	1
24"	1223011	1223211	1
30"	1323011	-	2
36"	1423011	-	2
42"	1823011	-	2
48"	1523011	-	2
513mm	1623011	1623211	1
1000mm	1723011	-	2

Width	Right Hinged	Left Hinged	# of Doors
12"	1028022	1028222	1
15"	1928022	1928222	1
18"	1128022	1128222	1
24"	1228022	1228222	1
30"	1328022	-	2
36"	1428022	-	2
42"	1828022	-	2
48"	1528022	-	2
513mm	1628022	1628222	1
1000mm	1728022	-	2

- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

End Sink Units



Description	54" Wide	58" Wide
Full Door Unit	9900019	9900020
False Panel Unit	9900023	9900024

- Full door end sink unit shown
- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

Sink Support Kit



Width	Item Number
12"	SSK1012
15"	SSK1015
18"	SSK1018
24"	SSK1024
30"	SSK1030
36"	SSK1036
42"	SSK1042
48"	SSK1048
54"	SSK1054
58"	SSK1058
Corner Unit	SSK0000

- Metal frame to support and position an undermount sink
- Sink Support Kit for Corner Unit is 26" deep and not exactly as shown

BASE CABINETS

31-3/4" HIGH X 22" DEEP

Open Units



Width	Item Number
12"	1080000
15"	1980000
18"	1180000
24"	1280000
30"	1380000
36"	1480000
42"	1880000
48"	1580000
513mm	1680000
1000mm	1780000

Full Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1080011	1080211	1
15"	1980011	1980211	1
18"	1180011	1180211	1
24"	1280011	1280211	1
30"	1380011	-	2
36"	1480011	-	2
42"	1880011	-	2
48"	1580011	-	2
513mm	1680011	1680211	1
1000mm	1780011	-	2

Door/Drawer Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1080022	1080222	1
15"	1980022	1980222	1
18"	1180022	1180222	1
24"	1280022	1280222	1
30"	1380022	-	2
36"	1480022	-	2
42"	1880022	-	2
48"	1580022	-	2
513mm	1680022	1680222	1
1000mm	1780022	-	2



Width	Item Number
30"	1380093
36"	1480093
42"	1880093
48"	1580093
1000mm	1780093



Width	Right Hand	Left Hand
30"	1380942	1380924
36"	1480942	1480924
42"	1880942	1880924
48"	1580942	1580924
1000mm	1780942	1780924

- Left hand door/drawer unit shown
- Right hand reversed



Width	Right Hand	Left Hand
30"	1380042	1380024
36"	1480042	1480024
42"	1880042	1880024
48"	1580042	1580024
1000mm	1780042	1780024

- Left hand door/drawer unit shown
- Right hand reversed

BASE CABINETS

31-3/4" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider

Width	Item Number
15"	1980090
18"	1180090
24"	1280090
30"	13800NN
36"	14800NN
42"	18800NN
48"	15800NN
513mm	1680090
1000mm	17800NN



up to 24" wide



30" and wider

Width	Item Number
12"	1080044
15"	1980044
18"	1180044
24"	1280044
30"	1380044
36"	1480044
42"	1880044
48"	1580044
513mm	1680044
1000mm	1780044



up to 24" wide



30" and wider

Width	Item Number
12"	1080055
15"	1980055
18"	1180055
24"	1280055
30"	1380055
36"	1480055
42"	1880055
48"	1580055
513mm	1680055
1000mm	1780055



Width	Item Number
30"	1380090
36"	1480090
42"	1880090
48"	1580090
1000mm	1780090



Width	Item Number
30"	1380097
36"	1480097
42"	1880097
48"	1580097
1000mm	1780097



Width	Item Number
30"	1380099
36"	1480099
42"	1880099
48"	1580099
1000mm	1780099

BASE CABINETS

31-3/4" HIGH X 22" DEEP

File Drawer Units



up to 24" wide



30" and wider

Width	Item Number
15"	1980100
18"	1180100
24"	1280100
30"	1380100
36"	1480100
42"	1880100
48"	1580100
513mm	1680100
1000mm	1780100



up to 24" wide



30" and wider

Width	Item Number
15"	1980111
18"	1180111
24"	1280111
30"	1380111
36"	1480111
42"	1880111
48"	1580111
513mm	1680111
1000mm	1780111

Sliding Door Units



Width	Item Number
30"	1380711
36"	1480711
42"	1880711
48"	1580711
1000mm	1780711

Waste/Storage Bin Units



up to 24" wide



30" and wider



Tilt Out Unit
(shown with
optional bin)



Pull Out Unit

Width	Tilt Out Units	Pull Out Units	Bin Capacity
15"	1980500	1980511	1
18"	1180500	1180511	1
24"	1280500	1280511	1
30"	1380500	1380511	2
36"	1480500	1480511	2
42"	1880500	1880511	2
48"	1580500	1580511	2
513mm	1680500	1680511	1
1000mm	1780500	1780511	2

- Add option code 84 to order stainless steel waste/storage bin, or add option code 85 for PVC waste/storage bin for Tilt Out units
- Tilt out units are hinged at the bottom
- Pull out units include full extension runners

BASE CABINETS

31-3/4" HIGH X 22" DEEP

Sink Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1083011	1083211	1
15"	1983011	1983211	1
18"	1183011	1183211	1
24"	1283011	1283211	1
30"	1383011	-	2
36"	1483011	-	2
42"	1883011	-	2
48"	1583011	-	2
513mm	1683011	1683211	1
1000mm	1783011	-	2



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
12"	1088022	1088222	1
15"	1988022	1988222	1
18"	1188022	1188222	1
24"	1288022	1288222	1
30"	1388022	-	2
36"	1488022	-	2
42"	1888022	-	2
48"	1588022	-	2
513mm	1688022	1688222	1
1000mm	1788022	-	2

- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coat (galvanized) inner horizontal bottom base plate as standard



End Sink Units



Description	54" Wide	58" Wide
Full Door Unit	9980017	9980018
False Panel Unit	9980021	9980022

- Full door end sink unit shown
- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

Pocket Door Sink Units

Base bottom in
down positionBase bottom in
knee space position


Width	Item Number
30"	138KS22
36"	148KS22
42"	188KS22
48"	158KS22
1000mm	178KS22

- Incorporates pocket doors and fold-up base to convert from sink cabinet design to knee space design
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet
- 42" wide cabinets and larger are ADA compliant

BASE CABINETS

31-3/4" HIGH X 22" DEEP


Diagonal Corner Sink Units



Description	Right Hinged	Left Hinged
Full Door Unit	1087211	1087511
False Panel Unit	1087722	1087622


- Left hinged, full door corner sink unit shown
- Add option code 90 to order Sink Hanger Supports
- Sink units are constructed with satin coated (galvanized) bottom, powder coated to match cabinet

Knee Space Sink Unit Dimensions



- Gable Legs are 1-1/2" wide
- Knee Space ships unassembled
- Accommodates a 6" deep sink
- 30" wide units with gable(s) are not ADA compliant
- See dotted line in figure to the right for ADA requirements

Knee Space Sink Unit



Width	Item Number
30"	13K3011
33"	10K3011
36"	14K3011
42"	18K3011
48"	15K3011
1000mm	17K3011

Center Panel Only



Width	Item Number
30"	13KP011
33"	10KP011
36"	14KP011
42"	18KP011
48"	15KP011
1000mm	17KP011

Left Leg With Center Panel



Width	Item Number
30"	13L3011
33"	10L3011
36"	14L3011
42"	18L3011
48"	15L3011
1000mm	17L3011

Right Leg With Center Panel



Width	Item Number
30"	13KR011
33"	10KR011
36"	14KR011
42"	18KR011
48"	15KR011
1000mm	17KR011

MOBILE BASE CABINETS

Mobile Cabinets



Our steel cabinets can be fitted with your choice of door and drawer fronts made from painted steel, stainless steel, wood, plastic laminate or phenolic resin in inset or architectural overlay design.

Heavy duty locking 3-1/2" twin wheel casters provide a total load rating (including mobile cabinet and work top) of 500lbs. Tip resisting counter weights are provided to aid in safe operation. Resulting mobility allows the user to quickly and easily move to another workstation, as tasks change. Use option code 77 for casters on cabinets.

OPTIONS

See Steel Specification & Options section for information on the following options:

- Door and drawer handles
- Door hinges
- Drawer body and suspension
- Lock and key options

BASE CABINET OPTIONS

Option	Description	Option	Description
00	Modified door/drawer configuration	74A	Toe kick - 1-1/2" on back
02	Label holders on all doors/drawers	75	Toe kick - front, back & two sides
03	Numeric or alpha identification disks on all doors/drawers (customer details required)	76	Replace toe kick with four square legs
04	Standard louver configuration on doors (6 louvers at 3" wide x 5-1/2" high)	77	Replace toe kick with four casters
05	Non-standard louver configuration on doors (customer drawing required)	78	Replace toe kick with customer specified casters
07	Cut-outs for electrical and/or plumbing services (customer drawing required)	79	Other customer specified base or toe kick
08	Security proof panel(s)	84	Stainless steel bin insert
09	Stainless steel waste flap	85	PVC bin insert
10	One fixed drawer divider per drawer (centered in drawer)	89	Lazy Susan storage unit for corner cabinets
11	Multiple drawer dividers per drawer (customer drawing required)	90	Sink supports - for under counter mounted sinks
12	Pull out writing board * - cold rolled steel	AD	Add - a - drawer design (6" high construction, applies to top drawer)
13	Pull out writing board * - cold rolled steel & high pressure laminate writing surface	BF	Replace toe kick with 4" bullet feet
14	Ventilation grill in toe kick base	C25	Replace toe kick with 250lb casters
15	Finished back on cabinet	C35	Replace toe kick with 350lb casters
16	Ship handles not installed	DC2	Magnetic door catch
17	Nylon levelers	F1	Welded finished top
17C	Nylon levelers with counter weight	F2	Welded finished top & back
18	Horizontally mounted handles on doors	F26	200lb full extension runners on file drawers
19	Full width drawer on 5" high apron	FCW	Front corners of doors and drawers welded
26	200lb full extension drawer runners	I01	100lb full extension runners with interlock
40	One extra cabinet shelf	I26	200lb full extension runners with interlock
41	Two extra cabinet shelves	LA	Bilingual label (French/English)
42	One interior cabinet shelf on full extension runners (replaces fixed shelf)	LF	Longitudinal file hangers on file drawers
43	Two interior cabinet shelves on full extension runners (replaces fixed shelves)	NB	No base, cabinet to sit on floor, levelers included
44	Three interior cabinet shelves on full extension runners (replaces fixed shelves)	NC	No base, cabinet to sit on 4" high curb, levelers not included
48	Keyboard tray (replaces drawer)	NF	No cabinet floor and toe kick is integral with door
50	304 stainless steel drawer bodies in power coated cold rolled steel cabinet	OL	Overlay design
51	Component made from 304 stainless steel	PC	Foot pedal cutouts centered
52	Component made from 316 stainless steel	PDC	Positive door catch (standard on overlay cabinets)
53	Galvenneal steel cabinet construction (Sigma Barrier™)	PFE	Pull out writing board * on full extension runners
54	304 stainless steel drawer bodies in galvenneal steel cabinet (Sigma Barrier™)	PL	Foot pedal cutouts left hand side
56	Other customer specified material	PR	Foot pedal cutouts right hand side
59	Customer specified runners (customer details required)	PW	Pre-wired (customer details required)
70	Front toe kick	RB	Removable back on drawer cabinet
71	Toe kick - front & left side	SK	Skid item
72	Toe kick - front & right side	S01	100lb full extension soft close drawer runners
73	Toe kick - front & 2 sides	S26	200lb full extension soft close drawer runners
74	Toe kick - front & back	T1	Replace one shelf with one pull out tray on full extension runners with 2" high lip
		T2	Replace two shelves with two pull out trays on full extension runners with 2" high lip
		T3	Replace three shelves with three pull out trays on full extension runners with 2" high lip
		WF	Wood fronts - see specifications
		WO	Wood fronts provided and installed by others

Note: * 30" wide cabinets and larger receive two pull out writing boards, one pull out board is provided on full width drawer cabinets.

All dimensions and sizes shown are nominal unless otherwise noted.
Specifications and details are based on product information at the time of printing and may change at any time without notice.
Mott Manufacturing reserves the right to change dimensions, specifications and manufacturing details at any time without notice.