

SUSPENDED CABINETS
24" HIGH X 22" DEEP
Open Units


Note: Add option code 15 to order a finished back on these units.

Width	Item Number
15"	1950000
18"	1150000
24"	1250000
30"	1350000
36"	1450000
42"	1850000
48"	1550000
513mm	1650000
1000mm	1750000

Full Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1950011	1950211	1
18"	1150011	1150211	1
24"	1250011	1250211	1
30"	1350011	-	2
36"	1450011	-	2
42"	1850011	-	2
48"	1550011	-	2
513mm	1650011	1650211	1
1000mm	1750011	-	2

- Cabinets receive one shelf

Door/Drawer Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1950022	1950222	1
18"	1150022	1150222	1
24"	1250022	1250222	1
30"	1350022	-	2
36"	1450022	-	2
42"	1850022	-	2
48"	1550022	-	2
513mm	1650022	1650222	1
1000mm	1750022	-	2

- Cabinets up to 24" wide receive one shelf, while cabinets 25" and wider do not have a shelf



Width	Right Hinged	Left Hinged
30"	1350041	1350014
36"	1450041	1450014
42"	1850041	1850014
48"	1550041	1550014
1000mm	1750041	1750014

- Left hand Door/Drawer unit shown, right hand reversed
- Cabinet receives one shelf

Sliding Door Units


Width	Item Number
30"	1350711
36"	1450711
42"	1850711
48"	1550711
1000mm	1750711

- Cabinets receive one shelf

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Drawer Units

Note: Add option code 15 to order a finished back on these units.



up to 24" wide



30" and wider

Width	Item Number
15"	1950055
18"	1150055
24"	1250055
30"	1350055
36"	1450055
42"	1850055
48"	1550055
513mm	1650055
1000mm	1750055



up to 24" wide



30" and wider

Width	Item Number
15"	1950066
18"	1150066
24"	1250066
30"	1350066
36"	1450066
42"	1850066
48"	1550066
513mm	1650066
1000mm	1750066



up to 24" wide



30" and wider

Width	Item Number
15"	1950044
18"	1150044
24"	1250044
30"	1350044
36"	1450044
42"	1850044
48"	1550044
513mm	1650044
1000mm	1750044



Width	Item Number
30"	1350097
36"	1450097
42"	1850097
48"	1550097
1000mm	1750097



Width	Item Number
30"	1350098
36"	1450098
42"	1850098
48"	1550098
1000mm	1750098

File Drawer Units



up to 24" wide



30" and wider

Width	Item Number
15"	1950122
18"	1150122
24"	1250122
30"	1350122
36"	1450122
42"	1850122
48"	1550122
513mm	1650122
1000mm	1750122

SUSPENDED CABINETS
24" HIGH X 22" DEEP
File Drawer Units

Note: Add option code 15 to order a finished back on these units.



up to 24" wide



30" and wider

Width	Item Number
15"	1950177
18"	1150177
24"	1250177
30"	1350177
36"	1450177
42"	1850177
48"	1550177
513mm	1650177
1000mm	1750177



up to 24" wide



30" and wider

Width	Item Number
15"	1950188
18"	1150188
24"	1250188
30"	1350188
36"	1450188
42"	1850188
48"	1550188
513mm	1650188
1000mm	1750188

Sink Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1953011	1953211	1
18"	1153011	1153211	1
24"	1253011	1253211	1
30"	1353011	-	2
36"	1453011	-	2
42"	1853011	-	2
48"	1553011	-	2
513mm	1653011	1653211	1
1000mm	1753011	-	2

Acid Storage Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1951011	1951211	1
18"	1151011	1151211	1
24"	1251011	1251211	1
30"	1351011	-	2
36"	1451011	-	2
42"	1851011	-	2
48"	1551011	-	2
513mm	1651011	1651211	1
1000mm	1751011	-	2

- Polyethylene liner and phenolic resin shelf are chemically resistant in acid cabinets
- Add option code VK to order a Corrugated Polyolefin Vent Kit for these units
- Note: Acid Cabinets must be vented. Mott recommends venting acid cabinets with a Polyolefin Vent Kit
- Contact Mott for more details

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Solvent Storage Units

Note: Add option code 15 to order a finished back on these units.



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1952211	1952311	1
18"	1152211	1152311	1
24"	1252211	1252311	1
30"	1352211	-	2
36"	1452211	-	2
42"	1852211	-	2
48"	1552211	-	2
513mm	1652211	1652311	1
1000mm	1752211	-	2

- Solvent units are not insulated
- Self closing doors are standard, use option FL to add a fusible link
- Solvent Units are designed to meet OSHA and NFPA 30 standards
- Contact Mott for more details

SUSPENDED CABINETS

18-3/4" HIGH X 22" DEEP

Open Units

Note: Add option code 15 to order a finished back on these units.



Width	Item Number
15"	1960000
18"	1160000
24"	1260000
30"	1360000
36"	1460000
42"	1860000
48"	1560000
513mm	1660000
1000mm	1760000

Full Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1960011	1960211	1
18"	1160011	1160211	1
24"	1260011	1260211	1
30"	1360011	-	2
36"	1460011	-	2
42"	1860011	-	2
48"	1560011	-	2
513mm	1660011	1660211	1
1000mm	1760011	-	2

Door/Drawer Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1960022	1960222	1
18"	1160022	1160222	1
24"	1260022	1260222	1
30"	1360022	-	2
36"	1460022	-	2
42"	1860022	-	2
48"	1560022	-	2
513mm	1660022	1660222	1
1000mm	1760022	-	2

SUSPENDED CABINETS
18-3/4" HIGH X 22" DEEP
Door/Drawer Units

Note: Add option code 15 to order a finished back on these units.



Width	Right Hinged	Left Hinged
30"	1360041	1360014
36"	1460041	1460014
42"	1860041	1860014
48"	1560041	1560014
1000mm	1760041	1760014

• Right hand Door/Drawer Unit shown, left hand reversed

Drawer Units


up to 24" wide



30" and wider

Width	Item Number
15"	1960044
18"	1160044
24"	1260044
30"	1360044
36"	1460044
42"	1860044
48"	1560044
513mm	1660044
1000mm	1760044

Drawer Units


up to 24" wide



30" and wider

Width	Item Number
15"	1960177
18"	1160177
24"	1260177
30"	1360177
36"	1460177
42"	1860177
48"	1560177
513mm	1660177
1000mm	1760177

Sink Units


up to 24" wide




30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1963011	1963211	1
18"	1163011	1163211	1
24"	1263011	1263211	1
30"	1363011	-	2
36"	1463011	-	2
42"	1863011	-	2
48"	1563011	-	2
513mm	1663011	1663211	1
1000mm	1763011	-	2

SUSPENDED COMPONENTS


Pencil Drawers and Pull Out Writing Tablets



Width	Pencil Drawer	Pull Out Writing Tablet
15"	PEN0015	PUL0015
18"	PEN0018	PUL0018
20"	PEN0020	PUL0020
24"	PEN0024	PUL0024
30"	PEN0030	PUL0030

- Pencil Drawers and Pull Out Writing Tablets suspends from tables with hanging rails and are supplied with full extension runners
- When ordering a table with aprons, add option code HU to suspend Pencil Drawer
- Pencil Drawers and Pull Out Writing Tablets do not fit on a 22" deep table

Hanging CPU Holder



Item Number

CPU1010

- Suspends from tables with hanging rails and is supplied with full extension runners to slide out CPU unit
- Inside clear dimension of CPU holder is 9-3/16"
- Hanging CPU Holder does not fit on a 22" deep table
- When ordering a table with aprons, add option code HU to suspend CPU holder

SUSPENDED CABINET OPTIONS

Option	Description
00	Modified door/drawer configuration
02	Label holders on all doors/drawers
03	Numeric or alpha identification disks on all doors/drawers (customer details required)
04	Standard louver configuration on doors (6 louvers at 3" wide x 5-1/2" high)
05	Non-standard louver configuration on doors (customer drawing required)
07	Cut-outs for electrical and/or plumbing services (customer drawing required)
08	Security proof panel(s)
09	Stainless steel waste flap
10	One fixed drawer divider per drawer (centered in drawer)
11	Multiple drawer dividers per drawer (customer drawing required)
12	Pull out writing board * - cold rolled steel
13	Pull out writing board - cold rolled steel & high pressure laminate writing surface
15	Finished back on cabinet
16	Ship handles not installed
18	Horizontally mounted handles on doors
26	200lb full extension drawer runners
40	One extra cabinet shelf
41	Two extra cabinet shelves
42	One interior cabinet shelf on full extension runners (replaces fixed shelf)
48	Keyboard tray (replaces drawer)
50	304 stainless steel drawer bodies in power coated cold rolled steel cabinet
51	Component made from 304 stainless steel
52	Component made from 316 stainless steel
53	Galvenneal steel cabinet construction (Sigma Barrier™)
54	304 stainless steel drawer bodies in galvenneal steel cabinet (Sigma Barrier™)
56	Other customer specified material
59	Customer specified runners (customer details required)

See Steel Specification & Options section for information on the following options:

- Door and drawer handles
- Door hinges
- Drawer body and suspension
- Lock and key options

Note: *30" wide cabinets and larger receive two pull out writing boards, one pull out board is provided on full width drawer cabinets.

Note: **Designed for use on solvent cabinets with self-closing. A fusible link holds doors open during use and melts at 165° F to automatically close the doors in the event of a fire.

Option	Description
84	Stainless steel bin insert
85	PVC bin insert
89	Lazy Susan storage unit for corner cabinets
90	Sink supports - for under counter mounted sinks
BL	BASE label (use this option code to replace standard label on acid cabinet)
CL	CORROSIVES label (use this option code to replace standard label on acid cabinet)
DC2	Magnetic door catch
F1	Welded finished top
F2	Welded finished top & back
F26	200lb full extension runners on file drawers
FCW	Front corners of doors and drawers welded
FL	Fusible link with self closer for insulated solvent unit**
I01	100lb full extension runners with interlock
I26	200lb full extension runners with interlock
LA	Bilingual label (French/English)
LF	Longitudinal file hangers on file drawers
OL	Overlay design (not available on Solvent Cabinets)
PP	Pre-plumbed (option code not to be used for fume hoods and customer details required)
PW	Pre-wired (customer details required)
PDC	Positive door catch (standard on overlay cabinets)
PFE	Pull out writing board * on full extension runners
RB	Removable back on drawer cabinet
SK	Skid item
SP	Silk screen printing (to replace standard labels on acid cabinet and/or flammable cabinet)
S01	100lb full extension soft close drawer runners
S26	200lb full extension soft close drawer runners
T1	Replace one shelf with one pull out tray on full extension runners with 2" high lip
T2	Replace two shelf with two pull out tray on full extension runners with 2" high lip
T3	Replace three shelf with three pull out tray on full extension runners with 2" high lip
VK	Corrugated Polyolefin vent kit for acid cabinets
VS	Vent holes on side of cabinet
WF	Wood fronts - see specifications
WO	Wood fronts provided and installed by others