

SUSPENDED CABINETS
24" HIGH X 22" DEEP
Open Units


Width	Item Number
15"	1950000W
18"	1150000W
24"	1250000W
30"	1350000W
36"	1450000W
42"	1850000W
48"	1550000W
513mm	1650000W
1000mm	1750000W

- Actual depths 7/8" less than nominal.

Full Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1950011W	1950211W	1
18"	1150011W	1150211W	1
24"	1250011W	1250211W	1
30"	1350011W	-	2
36"	1450011W	-	2
42"	1850011W	-	2
48"	1550011W	-	2
513mm	1650011W	1650211W	1
1000mm	1750011W	-	2

- Cabinets receive one shelf.

Door/Drawer Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1950022W	1950222W	1
18"	1150022W	1150222W	1
24"	1250022W	1250222W	1
30"	1350022W	-	2
36"	1450022W	-	2
42"	1850022W	-	2
48"	1550022W	-	2
513mm	1650022W	1650222W	1
1000mm	1750022W	-	2

- Cabinets up to 24" wide receive one shelf.
- Cabinets 25" and wider do not have a shelf.



Width	Right Hand	Left Hand
30"	1350041W	1350014W
36"	1450041W	1450014W
42"	1850041W	1850014W
48"	1550041W	1550014W
1000mm	1750041W	1750014 W

- Left hand Door/Drawer Unit shown.
- Right hand Door/Drawer Units are reverse of what is shown.
- Cabinets receive one shelf.

Sliding Door Units


Width	Item Number
30"	1350711W
36"	1450711W
42"	1850711W
48"	1550711W
1000mm	1750711W

- Actual depths 7/8" less than nominal.
- Cabinets receive one shelf.

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Drawer Units



up to 24" wide



30" and wider

Width	Item Number
15"	1950044W
18"	1150044W
24"	1250044W
30"	1350044W
36"	1450044W
42"	1850044W
48"	1550044W
513mm	1650044W
1000mm	1750044W



up to 24" wide



30" and wider

Width	Item Number
15"	1950055W
18"	1150055W
24"	1250055W
30"	1350055W
36"	1450055W
42"	1850055W
48"	1550055W
513mm	1650055W
1000mm	1750055W



up to 24" wide



30" and wider

Width	Item Number
15"	1950066W
18"	1150066W
24"	1250066W
30"	1350066W
36"	1450066W
42"	1850066W
48"	1550066W
513mm	1650066W
1000mm	1750066W



Width	Item Number
30"	1350097W
36"	1450097W
42"	1850097W
48"	1550097W
1000mm	1750097W



Width	Item Number
30"	1350098W
36"	1450098W
42"	1850098W
48"	1550098W
1000mm	1750098W

SUSPENDED CABINETS
24" HIGH X 22" DEEP
File Drawer Units


up to 24" wide



30" and wider



up to 24" wide



30" and wider



up to 24" wide



30" and wider

Sink Units


up to 24" wide



30" and wider

Acid Storage Units


up to 24" wide



30" and wider

Width	Item Number
15"	1950122W
18"	1150122W
24"	1250122W
30"	1350122W
36"	1450122W
42"	1850122W
48"	1550122W
513mm	1650122W
1000mm	1750122W

Width	Item Number
15"	1950177W
18"	1150177W
24"	1250177W
30"	1350177W
36"	1450177W
42"	1850177W
48"	1550177W
513mm	1650177W
1000mm	1750177W

Width	Item Number
15"	1950188W
18"	1150188W
24"	1250188W
30"	1350188W
36"	1450188W
42"	1850188W
48"	1550188W
513mm	1650188W
1000mm	1750188W

Width	Right Hinged	Left Hinged	# of Doors
15"	1953011W	1953211W	1
18"	1153011W	1153211W	1
24"	1253011W	1253211W	1
30"	1353011W	-	2
36"	1453011W	-	2
42"	1853011W	-	2
48"	1553011W	-	2
513mm	1653011W	1653211W	1
1000mm	1753011W	-	2

- Add option code 90 to order Sink Hanger Supports.

Width	Right Hinged	Left Hinged	# of Doors
15"	1951011W	1951211W	1
18"	1151011W	1151211W	1
24"	1251011W	1251211W	1
30"	1351011W	-	2
36"	1451011W	-	2
42"	1851011W	-	2
48"	1551011W	-	2
513mm	1651011W	1651211W	1
1000mm	1751011W	-	2

- Cabinet receives a molded liner and one chemically resistant shelf.
- Add option code VK to order a Corrugated Polyolefin Vent Kit for these units.
- Note: Acid Cabinets must be vented. Mott recommends venting acid cabinets with a Polyolefin Vent Kit.

SUSPENDED CABINETS

24" HIGH X 22" DEEP

Solvent Storage Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1952011W	1952211W	1
18"	1152011W	1152211W	1
24"	1252011W	1252211W	1
30"	1352011W	-	2
36"	1452011W	-	2
42"	1852011W	-	2
48"	1552011W	-	2
513mm	1652011W	1652211W	1
1000mm	1752011W	-	2

- Solvent Units are designed to meet OSHA and NFPA 30 standards.
- Self closing door(s) provided
- Add option code FL to order self closing door(s) with fusible link
- See page L10 for Solvent cabinet details.

SUSPENDED CABINETS

18" HIGH X 22" DEEP

Open Units



up to 24" wide

Width	Item Number
15"	1960000W
18"	1160000W
24"	1260000W
30"	1360000W
36"	1460000W
42"	1860000W
48"	1560000W
513mm	1660000W
1000mm	1760000W

- Actual depths 7/8" less than nominal.

Full Door Unit



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1960011W	1960211W	1
18"	1160011W	1160211W	1
24"	1260011W	1260211W	1
30"	1360011W	-	2
36"	1460011W	-	2
42"	1860011W	-	2
48"	1560011W	-	2
513mm	1660011W	1660211W	1
1000mm	1760011W	-	2

Door/Drawer Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1960022W	1960222W	1
18"	1160022W	1160222W	1
24"	1260022W	1260222W	1
30"	1360022W	-	2
36"	1460022W	-	2
42"	1860022W	-	2
48"	1560022W	-	2
513mm	1660022W	1660222W	1
1000mm	1760022W	-	2



Width	Right Hinged	Left Hinged
30"	1360041W	1360014W
36"	1460041W	1460014W
42"	1860041W	1860014W
48"	1560041W	1560014W
1000mm	1760041W	1760014W

- Right hand Door/Drawer Unit shown. Left hand reversed.

SUSPENDED CABINETS
18" HIGH X 22" DEEP
Drawer Units


up to 24" wide



30" and wider

Width	Item Number
15"	1960044W
18"	1160044W
24"	1260044W
30"	1360044W
36"	1460044W
42"	1860044W
48"	1560044W
513mm	1660044W
1000mm	1760044W

File Drawer Units


up to 24" wide



30" and wider

Width	Item Number
15"	1960177W
18"	1160177W
24"	1260177W
30"	1360177W
36"	1460177W
42"	1860177W
48"	1560177W
513mm	1660177W
1000mm	1760177W

Sink Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
15"	1963011W	1963211W	1
18"	1163011W	1163211W	1
24"	1263011W	1263211W	1
30"	1363011W	-	2
36"	1463011W	-	2
42"	1863011W	-	2
48"	1563011W	-	2
513mm	1663011W	1663211W	1
1000mm	1763011W	-	2

- Add option code 90 to order Sink Hanger Supports.

Hanging Rails


Width	Item Number
15"	HR00015
18"	HR00018
24"	HR00024
30"	HR00030
36"	HR00036
42"	HR00042
48"	HR00048
513mm	HR00020
1000mm	HR00040

- Pair of hanging rails only suitable for wood suspended cabinets.

SUSPENDED CABINET OPTIONS

Option Code	Description
00	Modified door/drawer configuration
02	Label holders on all doors/drawers
03	Numeric or alpha identification disks on all doors/drawers (customer details required)
04	Standard louver configuration on doors (6 louvers at 3" wide x 5-1/2" long)
05	Non-standard louver configuration on doors (customer drawing required)
07	Cut-outs for electrical and/or plumbing services (customer drawing required)
08	Security proof panel
09	Stainless steel waste flap
10	One fixed drawer divider per drawer (centered in drawer)
11	Multiple drawer dividers per drawer (customer drawing required)
12	Pull out writing board - all cold rolled steel*
13	Pull out writing board - high pressure laminate writing surface*
15	Finished back on cabinet
16	Ship handles not installed
18	Horizontally mounted handles
26	200 lbs full extension drawer runners
40	One extra cabinet shelf
41	Two extra cabinet shelves
42	One interior cabinet shelf on full extension runners (replaces fixed shelf)
43	Two interior cabinet shelves on full extension runners (replaces fixed shelf)
44	Three interior cabinet shelves on full extension runners (replaces fixed shelf)
48	Keyboard tray (replace drawer)
90	Sink supports - for under counter mounted sinks
BL	"BASE" label (use this option code to replace standard label on acid cabinets)
CL	CORROSIVES label
DC2	Magnetic Door Catch
F2	Finished top and back
FL	Fusible link with self closer for insulated unit**
LA	Bilingual label (French/English)
LF	Longitudinal file
OLD	Actual depth
PP	Pre-plumbed (option code not to be used for fume hoods)(customer details required)
PW	Pre-wired (option code not to be used for fume hoods) (customer details required)
PDC	Positive door catch
PFE	Pull out writing board on full extension runners*
RB	Removable back on drawer cabinet
SK	Skid item
SP	Silk screen printing (to replace standard labels on acid and/or flammable cabinets)
VK	Corrugated polyolefin vent kit for acid cabinets
VS	Vent holes on side of cabinet
T1	Replace one shelf with one pull out tray on full extension runners with 2" high lip
T2	Replace two shelves with two pull out trays on full extension runners with 2" high lip
T3	Replace three shelves with three pull out trays on full extension runners with 2" high lip

Notes: *30" wide cabinets and above receive 2 pull out writing boards unless it's a full width drawer cabinet.

Note: **Designed for use on solvent cabinets with self-closing and latching doors (option code 83). A fusible link holds doors open during use and melts at 165° F to automatically close doors in the event of a fire.