

WALL CABINETS
24" HIGH X 14" DEEP
Open Units


Width	Item Number
18"	5170010W
24"	5270010W
30"	5370010W
36"	5470010W
42"	5670010W
48"	5570010W
513mm	5970010W
1000mm	5870010W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5170030W	5171030W	1
24"	5270030W	5271030W	1
30"	5370030W	-	2
36"	5470030W	-	2
42"	5670030W	-	2
48"	5570030W	-	2
513mm	5970030W	5971030W	1
1000mm	5870030W	-	2

Hinged Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5170060W	5171060W	1
24"	5270060W	5271060W	1
30"	5370060W	-	2
36"	5470060W	-	2
42"	5670060W	-	2
48"	5570060W	-	2
513mm	5970060W	5971060W	1
1000mm	5870060W	-	2

Sliding Glass Door Units


Width	Item Number
30"	5370020W
36"	5470020W
42"	5670020W
48"	5570020W
1000mm	5870020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units


Width	Item Number
30"	5370040W
36"	5470040W
42"	5670040W
48"	5570040W
1000mm	5870040W

- Actual depth 7/8" less than nominal.

Sliding Door Units


Width	Item Number
30"	5370050W
36"	5470050W
42"	5670050W
48"	5570050W
1000mm	5870050W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Corner Units


Width	Door on Right	Door on Left
30"	5370330W	5370430W
36"	5470330W	5470430W
42"	5670330W	5670430W
48"	5570330W	5570430W
1000mm	5870330W	5870430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

WALL CABINETS

24" HIGH X 14" DEEP

Hinged Door Corner Units



Width	Door on Right	Door on Left
30"	5370360W	5370460W
36"	5470360W	5470460W
42"	5670360W	5670460W
48"	5570360W	5570460W
1000mm	5870360W	5870460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

WALL CABINETS

30" HIGH X 14" DEEP

Open Units



Width	Item Number
18"	5110010W
24"	5210010W
30"	5310010W
36"	5410010W
42"	5610010W
48"	5510010W
513mm	5910010W
1000mm	5810010 W

- Actual depth is 7/8" less than nominal.

Hinged Framed Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5110030W	5111030W	1
24"	5210030W	5211030W	1
30"	5310030W	-	2
36"	5410030W	-	2
42"	5610030W	-	2
48"	5510030W	-	2
513mm	5910030W	5911030W	1
1000mm	5810030W	-	2

Hinged Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5110060W	5111060W	1
24"	5210060W	5211060W	1
30"	5310060W	-	2
36"	5410060W	-	2
42"	5610060W	-	2
48"	5510060W	-	2
513mm	5910060W	5911060W	1
1000mm	5810060W	-	2

Sliding Glass Door Units



Width	Item Number
30"	5310020W
36"	5410020W
42"	5610020W
48"	5510020W
1000mm	5810020W

- Actual depth 7/8" less than nominal.

Sliding Framed Glass Door Units



Width	Item Number
30"	5310040W
36"	5410040W
42"	5610040W
48"	5510040W
1000mm	5810040W

- Actual depth 7/8" less than nominal.

WALL CABINETS
30" HIGH X 14" DEEP
Sliding Door Units


Width	Item Number
30"	5310050W
36"	5410050W
42"	5610050W
48"	5510050W
1000mm	5810050W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Corner Units


Width	Door on Right	Door on Left
30"	5310330W	5310430W
36"	5410330W	5410430W
42"	5610330W	5610430W
48"	5510330W	5510430W
1000mm	5810330W	5810430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

Hinged Door Corner Units


Width	Door on Right	Door on Left
30"	5310360W	5310460W
36"	5410360W	5410460W
42"	5610360W	5610460W
48"	5510360W	5510460W
1000mm	5810360W	5810460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

WALL CABINETS
36" HIGH X 14" DEEP
Open Units


Width	Item Number
18"	5120010W
24"	5220010W
30"	5320010W
36"	5420010W
42"	5620010W
48"	5520010W
513mm	5920010W
1000mm	5820010W

- Actual depth 7/8" less than nominal.

Hinged Framed Glass Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5120030W	5121030W	1
24"	5220030W	5221030W	1
30"	5320030W	-	2
36"	5420030W	-	2
42"	5620030W	-	2
48"	5520030W	-	2
513mm	5920030W	5921030W	1
1000mm	5820030W	-	2

Hinged Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5120060W	5121060W	1
24"	5220060W	5221060W	1
30"	5320060W	-	2
36"	5420060W	-	2
42"	5620060W	-	2
48"	5520060W	-	2
513mm	5920060W	5921060W	1
1000mm	5820060W	-	2

WALL CABINETS

36" HIGH X 14" DEEP

Sliding Glass Door Units



Width	Item Number
30"	5320020W
36"	5420020W
42"	5620020W
48"	5520020W
1000mm	5820020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units



Width	Item Number
30"	5320040W
36"	5420040W
42"	5620040W
48"	5520040W
1000mm	5820040W

- Actual depth 7/8" less than nominal.

Sliding Door Units



Width	Item Number
30"	5320050W
36"	5420050W
42"	5620050W
48"	5520050W
1000mm	5820050W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Corner Units



Width	Door on Right	Door on Left
30"	5320330W	5320430W
36"	5420330W	5420430W
42"	5620330W	5620430W
48"	5520330W	5520430W
1000mm	5820330W	5820430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

Hinged Door Corner Units



Width	Door on Right	Door on Left
30"	5320360W	5320460W
36"	5420360W	5420460W
42"	5620360W	5620460W
48"	5520360W	5520460W
1000mm	5820360W	5820460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

WALL CABINETS

48" HIGH X 14" DEEP

Open Units



Width	Item Number
18"	5130010W
24"	5230010W
30"	5330010W
36"	5430010W
42"	5630010W
48"	5530010W
513mm	5930010W
1000mm	5830010W

- Actual depth 7/8" less than nominal.

WALL CABINETS

48" HIGH X 14" DEEP

Hinged Framed Door Units



Width	Right Hinged	Left Hinged	# of Doors
18"	5130030W	5131030W	1
24"	5230030W	5231030W	1
30"	5330030W	-	2
36"	5430030W	-	2
42"	5630030W	-	2
48"	5530030W	-	2
513mm	5930030W	5931030W	1
1000mm	5830030W	-	2

up to 24" wide

30" and wider

Hinged Door Units



Width	Right Hinged	Left Hinged	# of Doors
18"	5130060W	5131060W	1
24"	5230060W	5231060W	1
30"	5330060W	-	2
36"	5430060W	-	2
42"	5630060W	-	2
48"	5530060W	-	2
513mm	5930060W	5931060W	1
1000mm	5830060W	-	2

up to 24" wide

30" and wider

Sliding Glass Door Units



Width	Item Number
30"	5330020W
36"	5430020W
42"	5630020W
48"	5530020W
1000mm	5830020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units



Width	Item Number
30"	5330040W
36"	5430040W
42"	5630040W
48"	5530040W
1000mm	5830040W

- Actual depth 7/8" less than nominal.

Sliding Door Units



Width	Item Number
30"	5330050W
36"	5430050W
42"	5630050W
48"	5530050W
1000mm	5830050W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Corner Units



Width	Door on Right	Door on Left
30"	5330330W	5330430W
36"	5430330W	5430430W
42"	5630330W	5630430W
48"	5530330W	5530430W
1000mm	5830330W	5830430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

WALL CABINETS

48" HIGH X 14" DEEP

Hinged Door Corner Units



Width	Door on Right	Door on Left
30"	5330360W	5330460W
36"	5430360W	5430460W
42"	5630360W	5630460W
48"	5530360W	5530460W
1000mm	5830360W	5830460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 14" wide on these units.

WALL CABINETS

24" HIGH X 18" DEEP

Open Units



Width	Item Number
18"	5180010W
24"	5280010W
30"	5380010W
36"	5480010W
42"	5680010W
48"	5580010W
513mm	5980010W
1000mm	5880010W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5180030W	5181030W	1
24"	5280030W	5281030W	1
30"	5380030W	-	2
36"	5480030W	-	2
42"	5680030W	-	2
48"	5580030W	-	2
513mm	5980030W	5981030W	1
1000mm	5880030W	-	2

Hinged Door Units



up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5180060W	5181060W	1
24"	5280060W	5281060W	1
30"	5380060W	-	2
36"	5480060W	-	2
42"	5680060W	-	2
48"	5580060W	-	2
513mm	5980060W	5981060W	1
1000mm	5880060W	-	2

Sliding Glass Door Units



Width	Item Number
30"	5380020W
36"	5480020W
42"	5680020W
48"	5580020W
1000mm	5880020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units



Width	Item Number
30"	5380040W
36"	5480040W
42"	5680040W
48"	5580040W
1000mm	5880040W

- Actual depth 7/8" less than nominal.

WALL CABINETS
24" HIGH X 18" DEEP
Sliding Door Units


Width	Item Number
30"	5380050W
36"	5480050W
42"	5680050W
48"	5580050W
1000mm	5880050W

- Actual depth is 7/8" less than nominal.

Hinged Framed Door Corner Units


Width	Door on Right	Door on Left
30"	5380330W	5380430W
36"	5480330W	5480430W
42"	5680330W	5680430W
48"	5580330W	5580430W
1000mm	5880330W	5880430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

Hinged Door Corner Units


Width	Door on Right	Door on Left
30"	5380360W	5380460W
36"	5480360W	5480460W
42"	5680360W	5680460W
48"	5580360W	5580460W
1000mm	5880360W	5880460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

WALL CABINETS
30" HIGH X 18" DEEP
Open Units


Width	Item Number
18"	5140010W
24"	5240010W
30"	5340010W
36"	5440010W
42"	5640010W
48"	5540010W
513mm	5940010W
1000mm	5840010W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5140030W	5141030W	1
24"	5240030W	5241030W	1
30"	5340030W	-	2
36"	5440030W	-	2
42"	5640030W	-	2
48"	5540030W	-	2
513mm	5940030W	5941030W	1
1000mm	5840030W	-	2

Hinged Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5140060W	5141060W	1
24"	5240060W	5241060W	1
30"	5340060W	-	2
36"	5440060W	-	2
42"	5640060W	-	2
48"	5540060W	-	2
513mm	5940060W	5941060W	1
1000mm	5840060W	-	2

WALL CABINETS

30" HIGH X 18" DEEP

Sliding Glass Door Units



Width	Item Number
30"	5340020W
36"	5440020W
42"	5640020W
48"	5540020W
1000mm	5840020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units



Width	Item Number
30"	5340040W
36"	5440040W
42"	5640040W
48"	5540040W
1000mm	5840040W

- Actual depth 7/8" less than nominal.

Sliding Door Units



Width	Item Number
30"	5340050W
36"	5440050W
42"	5640050W
48"	5540050W
1000mm	5840050W

- Actual depth 7/8" less than nominal.


Hinged Framed Door Corner Units



Width	Door on Right	Door on Left
30"	5340330W	5340430W
36"	5440330W	5440430W
42"	5640330W	5640430W
48"	5540330W	5540430W
1000mm	5840330W	5840430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

Hinged Door Corner Units



Width	Door on Right	Door on Left
30"	5340360W	5340460W
36"	5440360W	5440460W
42"	5640360W	5640460W
48"	5540360W	5540460W
1000mm	5840360W	5840460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

WALL CABINETS

36" HIGH X 18" DEEP

Open Units



Width	Item Number
18"	5150010W
24"	5250010W
30"	5350010W
36"	5450010W
42"	5650010W
48"	5550010W
513mm	5950010W
1000mm	5850010W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Units




Width	Right Hinged	Left Hinged	# of Doors
18"	5150030W	5151030W	1
24"	5250030W	5251030W	1
30"	5350030W	-	2
36"	5450030W	-	2
42"	5650030W	-	2
48"	5550030W	-	2
513mm	5950030W	5951030W	1
1000mm	5850030W	-	2

up to 24" wide

30" and wider

WALL CABINETS
36" HIGH X 18" DEEP
Hinged Door Units


up to 24" wide



30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5150060W	5151060W	1
24"	5250060W	5251060W	1
30"	5350060W	-	2
36"	5450060W	-	2
42"	5650060W	-	2
48"	5550060W	-	2
513mm	5950060W	5951060W	1
1000mm	5850060W	-	2

Sliding Glass Door Units


Width	Item Number
30"	5350020W
36"	5450020W
42"	5650020W
48"	5550020W
1000mm	5850020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units


Width	Item Number
30"	5350040W
36"	5450040W
42"	5650040W
48"	5550040W
1000mm	5850040W

- Actual depth 7/8" less than nominal.

Sliding Door Units


Width	Item Number
30"	5350050W
36"	5450050W
42"	5650050W
48"	5550050W
1000mm	5850050W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Corner Units


Width	Door on Right	Door on Left
30"	5350330W	5350430W
36"	5450330W	5450430W
42"	5650330W	5650430W
48"	5550330W	5550430W
1000mm	5850330W	5850430W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

Hinged Door Corner Units


Width	Door on Right	Door on Left
30"	5350360W	5350460W
36"	5450360W	5450460W
42"	5650360W	5650460W
48"	5550360W	5550460W
1000mm	5850360W	5850460W

- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

WALL CABINETS
48" HIGH X 18" DEEP
Open Units


Width	Item Number
18"	5160010W
24"	5260010W
30"	5360010W
36"	5460010W
42"	5660010W
48"	5560010W
513mm	5960010W
1000mm	5860010W

- Actual depth 7/8" less than nominal.

WALL CABINETS

48" HIGH X 18" DEEP

Hinged Framed Door Units




up to 24" wide

30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5160030W	5161030W	1
24"	5260030W	5261030W	1
30"	5360030W	-	2
36"	5460030W	-	2
42"	5660030W	-	2
48"	5560030W	-	2
513mm	5960030W	5961030W	1
1000mm	5860030W	-	2

Hinged Door Units



up to 24" wide

30" and wider

Width	Right Hinged	Left Hinged	# of Doors
18"	5160060W	5161060W	1
24"	5260060W	5261060W	1
30"	5360060W	-	2
36"	5460060W	-	2
42"	5660060W	-	2
48"	5560060W	-	2
513mm	5960060W	5961060W	1
1000mm	5860060W	-	2

Sliding Glass Door Units



Width	Item Number
30"	5360020W
36"	5460020W
42"	5660020W
48"	5560020W
1000mm	5860020W

- Actual depth 7/8" less than nominal.

Sliding Framed Door Units



Width	Item Number
30"	5360040W
36"	5460040W
42"	5660040W
48"	5560040W
1000mm	5860040W

- Actual depth 7/8" less than nominal.

Sliding Door Units



Width	Item Number
30"	5360050W
36"	5460050W
42"	5660050W
48"	5560050W
1000mm	5860050W

- Actual depth 7/8" less than nominal.

Hinged Framed Door Corner Units



Width	Door on Right	Door on Left
30"	5360330W	5360430W
36"	5460330W	5460430W
42"	5660330W	5660430W
48"	5560330W	5560430W
1000mm	5860330W	5860430W


- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

Hinged Door Corner Units




Width	Door on Right	Door on Left
30"	5360360W	5360460W
36"	5460360W	5460460W
42"	5660360W	5660460W
48"	5560360W	5560460W
1000mm	5860360W	5860460W


- Right door corner unit shown. Left door reversed.
- Blind portion of cabinet is 18" wide on these units.

WALL CABINETS
14"/18" DEEP
Hinged Framed Door Corner Units


Height	Right Hinged	Left Hinged	Depth
24"	5070030W	5071030W	14"
30"	5010030W	5011030W	14"
36"	5020030W	5021030W	14"
48"	5030030W	5031030W	14"
24"	5080030W	5081030W	18"
30"	5040030W	5041030W	18"
36"	5050030W	5051030W	18"
48"	5060030W	5061030W	18"

Hinged Door Corner Units


Height	Right Hinged	Left Hinged	Depth
24"	5070060W	5071060W	14"
30"	5010060W	5011060W	14"
36"	5020060W	5021060W	14"
48"	5030060W	5031060W	14"
24"	5080060W	5081060W	18"
30"	5040060W	5041060W	18"
36"	5050060W	5051060W	18"
48"	5060060W	5061060W	18"

Wall Cabinet Hanger Rails


Telescoping Wall Cabinet Hanger Rail



Fixed Width Wall Cabinet Hanger Rail

Width	Sigma Flex™		Sigma Frame™	
	Telescoping Rail	Fixed Width Rail	Telescoping Rail	Fixed Width Rail
24"	TWR0024	FWR0024	AWR0024	SWR0024
30"	TWR0030	FWR0030	AWR0030	SWR0030
36"	TWR0036	FWR0036	AWR0036	SWR0036
42"	TWR0042	FWR0042	AWR0042	SWR0042
48"	TWR0048	FWR0048	AWR0048	SWR0048
54"	TWR0054	FWR0054	AWR0054	SWR0054
60"	TWR0060	FWR0060	AWR0060	SWR0060
72"	TWR0072	FWR0072	AWR0072	SWR0072

- Wall Cabinet Hanger Rails are required to hang above-counter storage cabinets.
- Order rails to match width of support structure frame (Telescoping Rails are not suitable for cabinets smaller than 24" wide).
- Telescoping Rail is adjustable up to 5".
- Wall cabinets must be ordered with "J" Mounting Rails to hang on hanger rails.

Wall Cabinet "J" Mounting Rails


Width	Item Number
24"	JR00024
30"	JR00030
36"	JR00036
42"	JR00042
48"	JR00048
54"	JR00054
60"	JR00060
66"	JR00066
72"	JR00072

- A "J" Mounting Rail is required to hang metal above-counter storage cabinets on to Sigma Flex™ or Sigma Frame™ products.
- Order "J" rails to match width of support structure frame and Wall Cabinet Hanger Rail.
- "J" Rail ships loose.

WALL CABINET OPTIONS

Option Code	Description
00	Modified door/drawer configuration
02	Label holders on all doors/drawers
03	Numeric or alpha identification disks on all doors/drawers (customer details required)
09	Stainless steel waste flap
16	Ship handles not installed
18	Horizontally mounted handles
40	One extra cabinet shelf
41	Two extra cabinet shelves
45	DuoLite safety glass
46	Tempered glass
47	Customer specified glass
4A	Glass Tint 44% (Grey Smoke)
DS	Double sided wall cabinet
DC2	Magnetic door catch
OLD	Actual depth
SK	Skid item
SL	Seismic lip