

Adaptable Tables



1 CAMBRIDGE SERIES TABLES

tables

48" wide x 30" deep
60" wide x 30" deep
72" wide x 30" deep

- 1" black epoxy top included; 2" side/front overhangs, flush with the back of the table frame
- table leg opening (subtract 8.5" from top width)
- can accommodate lower shelf

2 CAMBRIDGE SHELVING KITS - STRUT

overhead strut shelving kit

48" table shelving kit
60" table shelving kit
72" table kit shelving

shelving kit includes:

- (2) 11 ga. painted steel strut uprights
- (2) shelf support brackets
- (2) bookstop rails
- (2) black plastic laminate shelves
- (2) bookend brackets
- bolts and hardware
- please note: electrical strip sold separately

3 ELECTRICAL STRIP ASSEMBLIES - STRUT

power strip and rail

48" electrical strip kit
60" electrical strip kit
72" electrical strip kit

- raceway with 15' grounded cord, on/off switch, and 15A 120V outlets mounted to a steel support rail
- 36" power strip for 48" tables; 12 outlets
- 48" power strip for 60" tables; 16 outlets
- 60" power strip for 72" tables; 20 outlets



tables

construction: fully welded frame with stainless steel bolted leg inserts

adjustable with 1" work surface

- 25"h to 37"h on levelers
- 30"h to 42"h on casters

mobility: leveling feet and casters are included if desired and specified

upper flexibility: optional shelving system bolts to rear table legs (adding 1.625" d)

lower flexibility: hanging rails support under work surface components (see pg. 9 & 10)

weight capacity: 1,000 lb. evenly distributed on frame (please account for top, storage & bench top equipment)

frame color: bright white or warm gray

strut shelving kits

adjustment: brackets are adjustable in height along the strut (typically rising 47" above the top)

shelving: 1" black plastic laminate banded with a durable 3mm black PVC edge

*phenolic shelving available upon request (see pg. 14)

flexibility: add/remove/modify kit on any standard Cambridge Series base table

weight capacity: supports 200 lb. evenly distributed per shelf

steel frame components: bright white or warm gray

shelving: black or gray

CAMBRIDGE SHELVEING KIT - SLOTTED SINGLE



single bay overhead slotted shelving

- 48" table shelving kit
- 60" table shelving kit
- 72" table shelving kit

shelving kit includes:

- (1) left slotted upright with service holes
- (1) right slotted upright
- (1) cross brace
- (2) painted steel shelves w/rear lip
- (4) tapered shelf brackets adjustable in 1" increments
- color: bright white only

BW

SERVICE VALVES



needle valve (left)

- VAC (yellow) - 61997762028
- AIR (orange) - 61997882021
- NIT (brown) - 61997812023

ADA ball valve (right)

- VAC (yellow) - 61997782028
- AIR (orange) - 61997872021
- NIT (brown) - 61997822023

BW

POST PLUMBING KIT



plumbing kit

54" post plumbing kit

plumbing kit includes:

- (3) color-coded quick release plugs at top of hose
- soft PEX braided stainless steel hoses for AIR/NIT with color coded PVC labels
- Teflon braided stainless steel hose for VAC with color coded PVC label
- (3) bright white blind caps

ELECTRICAL STRIP ASSEMBLIES - SLOTTED



economy

- 48" electrical strip kit
- 60" electrical strip kit
- 72" electrical strip kit

standard

- 48" electrical strip kit
- 60" electrical strip kit
- 72" electrical strip kit

premium

- 48" electrical strip kit
- 60" electrical strip kit
- 72" electrical strip kit

- twist lock plug on a 10' white cord and 90° cord grip
- raceway can be flipped so cord can be used either side of bench
- anodized aluminum with removable cover
- economy: (5) round standard 20A duplex outlets
- standard: (5) rectangular designer 20A duplex outlets, (2) type B USB ports
- premium: (6-8) rectangular designer 20A duplex outlets, (2) type B USB ports
- color: bright white only

BW

QUICK CONNECTS & HOSES



quick connects (left)

- VAC (yellow) - 25980222028
- AIR (orange) - 25980222021
- NIT (brown) - 25980222023

white hoses (right)

- VAC (yellow) - 25980402028
- AIR (orange) - 25980402021
- NIT (brown) - 25980402023

MONITOR ARM



monitor arm

IN-7000-FM

- supplied with mounting hardware for 6 different mounting methods
- color: silver, additional colors available



FLAT WHITE
248



SILVER
124



VISTA BLACK
104



desk edge
mount



through-hole
mount



wall mount



reverse wall/table
upright mount



side bolt
mount



grommet hole
mount



monitor arm

IN-9112-S-FM

CONFORM MONITOR ARM



single monitor arm

WR-1SDA-CCG



dual monitor arms

WR-2SDA-TPCCG

- single arm supports 1 computer monitor
- single arm maximum monitor height: 16"
- dual arm supports 2 computer monitors
- dual arm maximum monitor height: 19.5"
- color: silver
- maintenance free gas spring
- supports 5-17 lb. monitors
- integrated wire managers
- ability to fold down
- articulation range: 13"
- tilt range: -45° +90°
- clamp base

KEYBOARD TRAY



articulating slide-out tray

KT12-27

- size: 27"w
- 12.5" height adjustment range (7" up, 5.5" down)
- 23" track included
- lift to desired position and release

KEYBOARD & MOUSE TRAY



keyboard tray w/swivel mouse tray

KT12-19

- size: 19"w
- 12.5" height adjustment range (7" up, 5.5" down)
- 23" track included
- lift to desired position and release



TASK LIGHT



task light

- 12-5/16" wide light, 26 LEDs
- 21-5/16" wide light, 45 LEDs
- 22-9/16" wide light, 51 LEDs
- 34-3/8" wide light, 76 LEDs
- 46-3/16" wide light, 100 LEDs
- 58" wide light, 126 LEDs

- color: white
- screw mounting clips, in-line on/off rocker switch and power cord included
- 120 volts
- 4000K standard light temperature
- high efficiency with 87 Lm/W
- accessories sold separately

accessories

- 2) connection cord
- 3) direct connector
- 4) magnetic mounting clips



CPU HOLDER SUSPENDED



cpu holder

8-1/4" wide - CPU29157

- color: black
- adjustable from 8-1/4"w - 3-3/4"w and 20-1/2"h - 12"h

ANTI-FATIGUE MAT



floor mat

WR-AFM-2918

- color: black
- 29.5"w x 18"d x 3/4"h
- 100% polyurethane foam with puncture resistant surface
- patent pending foot hold for convenient repositioning
- 18" beveled edges



NOVUS TOOLBAR



table size	width	surface attachment
4'	40"	clamp mount*
5'	50"	clamp mount*
6'	63"	clamp mount*

- flexible organizer for your desktop with a variety of accessories to support your workspace
- (1) aluminum slatwall element with two slots for attaching accessories, monitor arms and task lights
- (2) aluminum posts with clamp mount



*available with cable port mount

NOVUS UPRIGHT BAR KIT




table size	width	attachment
4'	40.12"	mounting brackets
5'	52.13"	mounting brackets
6'	64.13"	mounting brackets

- (1) extruded aluminum slat rail bar with (2) mounting brackets and (2) end caps
- flush mount: attach bar to the back of the brackets - end caps not needed
- forward mount: attach bar to the face of the bracket and cap the ends

***For use with the Cambridge Table slotted upright system. Uprights not included.**

ACCESSORIES



my arm monitor arm
NV-MA-CM



LifTEC arm III
NV-LT-3



slatwall monitor arm
NV-CLII-209



multi storage board
NV-MSB-99



cell phone holder
NV-PHN-79



tablet holder
NV-TH-89



magnetic strip
NV-MS-109



pen holder
NV-PH-69



memo holder
NV-MB-39



pura-line pen cup
NV-PBH-69



penda small multi-purpose box
NV-SB-49



shelf
NV-SWSHV-600 23.62"w x 12"d
NV-SWSHV-600 31.5" w x 12"d
NV-SWSHV-600 39.37"w x 12"d

Painted Steel Cabinets



1 3-HIGH DRAWER SUSPENDED

full width drawers

15" wide - 1950044

18" wide - 1150044

24" wide - 1250044

- may interfere with caster swivel at lower table height adjustments

2 HINGED GLASS WALL CABINETS

double door

30" wide - 5310030

36" wide - 5410030

- 1/4" clear float glass
- (2) adjustable steel shelves included
- mounting channel included

3 BASE DOOR/DRAWER CABINETS

double door

30" wide - 1310022

36" wide - 1410022

48" wide - 1510022

- (1) adjustable steel shelf included
- standing height (34-3/4"h)

suspended cabinets

size: cabinets are 24"h and 22"d

construction: inset steel

drawers slides: full extension 100 lb. soft close slides

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges

finished back and tops

color: bright white or warm gray

wall cabinets

size: cabinets are 30"h and 14"d

construction: inset steel

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges

color: bright white or warm gray



DOOR/DRAWER SUSPENDED



left hinged (shown)

18" wide - 1150222

24" wide - 1240222

right hinged

18" wide - 1150022

24" wide - 1250022

- may interfere with caster swivel at lower table height adjustments

3-HIGH DRAWER SUSPENDED



full width drawers

15" wide - 1950044

18" wide - 1150044

24" wide - 1250044

- may interfere with caster swivel at lower table height adjustments

HINGED WALL CABINETS



steel double door

30" wide - 5310060

36" wide - 5410060



glass double door

30" wide - 5310030

36" wide - 5410030

- 1/4" clear float glass (glass wall cabinets only)
- (2) adjustable steel shelves included

DETERMINING DIMENSIONS



wall and suspended cabinets

size:

wall cabinets: 30"h x 14"d

suspended cabinets: 24"h x 22"d

construction: inset steel

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges

color: bright white or warm gray

suspended cabinets:

drawer slides: full extension 100 lb. softclose slides

hanging (suspended file cabinets): files hang front to back in the drawer

mounting channels included on wall shelves