



1 EPOXY PEGBOARD

sizes

- 24" wide x 30" high - 39 pegs - EPB2430
- 30" wide x 30" high - 46 pegs - EPB3030
- 30" wide x 36" high - 56 pegs - EPB3036

- material: 1" thick black (or gray) epoxy resin with beveled edges
- installation: 4 screw holes provided for wall attachment
- drip trough: 2" d stainless steel with 36" clear flexible drain tube
- some assembly required

2 MIXING FAUCET

mixing faucet (with VB)

L414-8VB-55-BO140S

- mixing faucet: deck-mounted, with vacuum breaker for hot and cold water
- handles: forged brass four-arm handles with color-coded index discs
- finish: forged brass with polished chrome plated finish

faucet

valve units: self-contained compression valve with replaceable stainless steel seats

gooseneck: 8" rigid/swing convertible gooseneck

inlet: 1-3/16" male shank with (2) 3/8" OD flexible copper tubes; furnished complete with locknut and washer

outlet: 3/8" NPS female outlet; included ten serration hose end or aerator

vacuum breaker (VB): integral vacuum breaker to prevent backflow; breaker has replaceable seat and ultra-light float cup to resist spilling at low flow

3 BLACK EPOXY RESIN

work surfaces

- 30" x 36" - EPX3036
- 30" x 48" - EPX3048
- 30" x 60" - EPX3060
- 30" x 72" - EPX3072
- 30" x 96" - EPX3096
- 30" x 100" - EPX30100
- 36" x 100" - EPX36100

sink top*

- 30" x 50" - EPX3050E
- 30" x 60" - EPX3060E
- 30" x 72" - EPX3072E

backsplash

- 4" x 60" - EPX0460
- 4" x 100" - EPX04100

sidesplash

- 4" x 29" - EPX0429

work surface

1" thick monolithic epoxy resin

edges: 1/8" bevel with drip groove on underside


splash: applied

adhesive: two-part epoxy provided as needed

other sizes can be readily fabricated; please consult your representative for lead times

*sinktop E: E indicates that the top includes a black D55 dop-in epoxy sink (25"x15"x10"), cutout and outlet; sinks are centred right to left in the work surface

STAINLESS PEGBOARDS



sizes

18" wide x 24" high - 15 pegs - V1824
 24" wide x 30" high - 32 pegs - V2430
 30" wide x 36" high - 40 pegs - V3036

- material: 20 gauge, type 304 stainless steel, polished to #4 finish
- installation: wall mount bracket shipped with each pegboard
- 6" long 1/2" diameter white polypropylene pegs
- cleaning: safe to autoclave
- integrated drip tray

see page 20 for epoxy pegboards

EPOXY SINK



sink

D55

- drop-in
- 25" w x 15" h x 10" d
- cutout and outlet
- colors: black and gray



outlet

OUTL-S03R-001

sink strainer

STRN-01-001

- for use with D55 sink
- colors: black and gray

STAINLESS STEEL SINK



sink

SLXD-2225

- 22" w x 16" h x 12" d
- cutout and strainer

ACRYLIC SPLASHGUARD



sizes

30" x 12" x 0.5" - ASGN3012-SS29-SS8
 30" x 18" x 0.75" - ASGN3018-SS29-SS14

- notched with stainless steel U-channel
- eased and flame-polished edges
- 2" radius upper corner

