

Painted Steel Cabinets



1 3-HIGH DRAWER SUSPENDED

full width drawers

15" wide - 1950044

18" wide - 1150044

24" wide - 1250044

Please note: may interfere with caster swivel at lower table height adjustments

BW

2 HINGED GLASS WALL CABINETS

double door

30" wide - 5310030

36" wide - 5410030

- 1/4" clear float glass
- (2) adjustable steel shelves included
- mounting channel included

BW

3 BASE DOOR/DRAWER CABINETS

double door

30" wide - 1310022

36" wide - 1410022

48" wide - 1510022

- (1) adjustable steel shelf included
- mounting channel included
- standing height (34-3/4"h)

BW

suspended cabinets

size: cabinets are 24"h and 22"d

construction: inset steel

drawers slides: 3/4 extension 100 lb. nylon roller slides

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges

finished back and tops

color: warm gray or bright white

wall cabinets

size: cabinets are 30"h and 14"d

construction: inset steel

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges

color: warm gray or bright white

mott
MANUFACTURING
ISO 9001:2015 REGISTERED

HINGED WALL CABINETS



left hinged (shown)

18" wide - 5111060

right hinged

18" wide - 5110060

- right hinged is reversed of what is shown
- (2) adjustable steel shelves included

BW

HINGED WALL CABINETS



double door

30" wide - 5310060

36" wide - 5410060

- (2) adjustable steel shelves included

BW

HINGED GLASS WALL CABINETS



left hinged (shown)

18" wide - 5111030

right hinged

18" wide - 5110030

- 1/4" clear float glass
- (2) adjustable steel shelves included

BW

DOOR/DRAWER SUSPENDED



left hinged (shown)

18" wide - 1150222

24" wide - 1240222

right hinged

18" wide - 1150022

24" wide - 1250022

- may interfere with caster swivel at lower table height adjustments

BW

FILE DRAWER SUSPENDED



box/file drawers

18" wide - 1150177

- may interfere with caster swivel at lower table height

wall and suspended cabinets

size:

wall cabinets: 30"h x 14"d

suspended cabinets: 24"h x 22"d

construction: inset steel

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges

color: warm gray, limited pieces available in bright white

suspended cabinets:

drawer slides (suspended cabinets): 3/4 extension 100 lb. nylon roller slides

hanging (suspended file cabinets): files hang front to back in the drawer

mounting channels included on wall shelves

see page 9 for additional suspended and wall cabinets



FULL DOOR CABINETS



double door

36" wide - 1410011

- (1) adjustable steel shelf included

BW

DOOR/DRAWER CABINETS



left hinged (shown)

18" wide - 1110222

24" wide - 1210222

right hinged

18" wide - 1110022

24" wide - 1210022

- right hinged door/drawer is reversed of what is shown

BW

SINK CABINETS



double door

36" wide - 1418022

48" wide - 1518022

- sink supports available for use with undermount sinks
- powdercoat galvanized steel on cabinet bottom

BW

END SINK CABINETS



double door

58" wide - 9900022

- sink supports available for use with undermount sinks
- powdercoat galvanized steel on cabinet bottom

BW

COMBINATION UNITS



left hinged (shown)

36" wide - 1410014

right hinged

36" wide - 1410041

- right hinged combination is reversed of what is shown
- (1) adjustable shelf included

FILE DRAWER SITTING HEIGHT



box/box/file drawers

18" wide - 1120188

- files hang front to back
- size: 28.75"h x 22"d

mott
MANUFACTURING
ISO 9001:2015 REGISTERED

3-HIGH DRAWER CABINETS



up to 24" wide



30" wide

full width drawers

18" wide - 1110090

24" wide - 1210090

30" wide - 1310090

BW

BW

BW

full width drawers - sitting height (28-3/4")

18" wide - 1120044-SH

24" wide - 1220044-SH

4-HIGH DRAWER CABINETS



up to 24" wide



30" wide and up

full width drawers

18" wide - 1110044

24" wide - 1210044

30" wide - 1310097

36" wide - 1410097

BW

DETERMINING DIMENSIONS



base cabinets

size: 34.75"h x 22"d

construction: inset steel

drawers slides: 3/4 extension 100 lb. nylon roller slides

hardware: aluminum wire pulls and stainless steel 5-knuckle hinges


leveling method: adjustable leveling feet

color: warm gray; most pieces also available in bright white

see page 9 for additional base cabinets



LEFT FRONT FILLER PANELS

**filler 34.75" high (left scribe)**

1" wide - FFL1001

4" wide - FFL1004

12" wide - FFL1012

filler 28.75" high (left scribe)

4" wide - FFL2004

toe base corner filler

9900001

filler panel

2" wide - FFC1002

6" wide - FFC1006

- toe base filler required for corners
- can be cut to meet wall conditions
- powder coated steel for inset cabinets
- color: warm gray and bright white

BW

RIGHT FRONT FILLER PANELS

**filler 34.75" high (right scribe)**

1" wide - FFR1001

4" wide - FFR1004

12" wide - FFR1012

filler 28.75" high (right scribe)

4" wide - FFR2004

- can be cut to meet wall conditions
- powder coated steel for inset cabinets
- color: warm gray and bright white

BW

REAR FILLER PANELS

**filler 34.75" high (rear scribe)**

1" wide - 24" deep tops - ECS1001

7" wide - 30" deep tops - ECS1007 BW

11" wide - hood cabinets - ECS1011 BW

14" wide - hood cabinets - ECS1014 BW

filler 34.75" high

14" wide - island ends - ECC1014 BW

filler 28.75" high (rear scribe)

1" wide - 24" deep tops - ECS2001

7" wide - 30" deep tops - ECS2007

- rear filler can be used for right/left hand applications
- scribe not required at island/peninsula ends
- can be cut to meet wall conditions
- powder coated steel for inset cabinets
- color: warm gray; limited pieces in bright white

WALL FILLER PANELS

**30" deep fillers**

1" wide - WFS6301

4" wide - WFS6304

15" wide - WFS6315

- for use with 14" d wall cabinets
- wall fillers can be used for right/left hand applications
- can be cut to meet wall conditions
- powder coated steel for inset cabinets
- color: warm gray and bright white

BW

KNEESPACE PANELS

**knee panel 34.75" high**

30" wide - SPC1030 BW

36" wide - SCP1036 BW

knee panel 28.75" high

30" wide - SPC2030

36" wide - SCP2036

- assembly includes single angles for attachment to adjacent cabinets
- upper panel is removable to allow access to the services
- lower panel stays fixed in order to apply the coved base
- overlapping lower and upper panel allows for adjustment on uneven floors
- powder coated steel for inset cabinets
- color: warm gray; limited pieces in bright white

FRONT RAILS

**3.75" high**

30" wide - FRA0030 BW

36" wide - FRA0036 BW

2" high

30" wide - FRB0030

- front rails can be used at kneespace locations
- powder coated steel for inset cabinets
- color: warm gray; limited pieces in bright white