

Adaptable Tables



1 CAMBRIDGE SERIES TABLES

tables

- 36" wide x 30" deep
- 48" wide x 30" deep
- 60" wide x 30" deep
- 72" wide x 30" deep
- *96" wide x 30" deep

- 1" black epoxy top included; 2" side/front overhangs, flush with the back of the table frame
- table leg opening (subtract 8.5" from top width)
- can accommodate lower shelf

BW

tables

construction: fully welded frame with stainless steel bolted leg inserts

adjustment: adjustable with 1" work surface

- 25"h to 37"h on levelers
- 30"h to 42"h on casters

mobility: leveling feet and casters are included if desired and specified

upper flexibility: optional shelving system bolts to rear table legs (adding 1.625"d)

lower flexibility: hanging rails support under work surface components (see pg. 9 & 10)

weight capacity: 1,000 lb. evenly distributed on frame (please account for top, storage & bench top equipment)

* 96" table supports 500 lb. evenly distributed

frame color: warm gray or bright white

2 CAMBRIDGE SHELVING KITS - STRUT

overhead strut shelving kit

- 48" table shelving kit
- 60" table shelving kit
- 72" table kit shelving
- 96" table kit shelving

shelving kit includes:

- (2) 11 ga. painted steel strut uprights
- (2) shelf support brackets
- (2) bookstop rails
- (2) black plastic laminate shelves
- (2) bookend brackets
- bolts and hardware
- please note: electrical strip sold separately

BW

strut shelving kits

adjustment: brackets are adjustable in height along the strut (typically rising 47" above the top)

shelving: 1" black plastic laminate banded with a durable 3mm black PVC edge
*phenolic shelving available upon request (see pg. 14)

flexibility: add/remove/modify kit on any standard Cambridge Series base table

weight capacity: supports 200 lb. evenly distributed per shelf

steel frame components: warm gray or bright white

shelving: black or gray

3 ELECTRICAL STRIP ASSEMBLIES - STRUT

power strip and rail

- 48" electrical strip kit
- 60" electrical strip kit
- 72" electrical strip kit
- 96" electrical strip kit

- raceway with 15' grounded cord, on/off switch, and 15A 120V outlets mounted to a steel support rail
- 36" power strip for 48" tables; 12 outlets
- 48" power strip for 60" tables; 16 outlets
- 60" power strip for 72" tables; 20 outlets
- 72" power strip for 96" tables; 24 outlets



BW



CAMBRIDGE SHELVEING KIT - SLOTTED SINGLE



single bay overhead slotted shelving

48" table shelving kit

60" table shelving kit

72" table shelving kit

shelving kit includes:

- (1) left slotted upright with service holes
- (1) right slotted upright
- (1) cross brace
- (2) painted steel shelves w/rear lip
- (4) tapered shelf brackets adjustable in 1" increments

BW

CAMBRIDGE SHELVEING KIT - SLOTTED DOUBLE



double bay overhead slotted shelving

60" table shelving kit

72" table shelving kit

shelving kit includes:

- (1) left slotted upright with service holes
- (1) right slotted upright
- (2) cross braces
- (4) painted steel shelves w/rear lip
- (8) tapered shelf brackets adjustable in 1" increments

BW

POST PLUMBING KIT



plumbing kit

54" post plumbing kit

plumbing kit includes:

- (3) color-coded quick release plugs at top of hose
- soft PEX braided stainless steel hoses for AIR/NIT with color coded PVC labels
- Teflon braided stainless steel hose for VAC with color coded PVC label
- (3) Bright White blind caps

SERVICE VALVES



needle valve (left)

VAC (yellow) - 61997762028

AIR (orange) - 61997882021

NIT (brown) - 61997812023

ADA ball valve (right)

VAC (yellow) - 61997782028

AIR (orange) - 61997872021

NIT (brown) - 61997822023

BW

QUICK CONNECTS & HOSES



quick connects (left)

VAC (yellow) - 25980222028

AIR (orange) - 25980222021

NIT (brown) - 25980222023

white hoses (right)

VAC (yellow) - 25980402028

AIR (orange) - 25980402021

NIT (brown) - 25980402023

ELECTRICAL STRIP ASSEMBLIES - SLOTTED



economy

48" electrical strip kit

60" electrical strip kit

72" electrical strip kit

standard

48" electrical strip kit

60" electrical strip kit

72" electrical strip kit

premium

48" electrical strip kit

60" electrical strip kit

72" electrical strip kit

- twist lock plug on a 10' white cord and 90° cord grip
- raceway can be flipped so cord can be used either side of bench
- anodized aluminum with removable cover
- economy: (5) round standard 20A duplex outlets
- standard: (5) rectangular designer 20A duplex outlets, (2) type B USB ports
- premium: (6-8) rectangular designer 20A duplex outlets, (2) type B USB ports

KEYBOARD TRAY



articulating slide-out tray*
IN-8491-8495

- size: 27" w
- 12.5" height adjustment range (7" up, 5.5" down)
- 23" track included
- lift to desired position and release

STANDARD KEYBOARD TRAY



keyboard tray
WR-ST185

ultra-thin keyboard tray
WR-ST185-UT

- mouse: right/left reversible with retainer
- smooth nylon glide track

BANANA BOARD



keyboard tray
WR-UB2180

- non-skid keyboard pads
- built-in document holder
- integrated cable management
- sliding reversible mouse surface
- mouse-trap mouse retainer
- maintenance-free ball bearing track

REVO KEYBOARD TRAY



keyboard tray w/mouse pad
WR-UB2100

- 0.25" solid phenolic construction
- built-in document holder
- integrated cable management
- mouse: right/left reversible with retainer
- maintenance-free ball bearing track

PINNACLE KEYBOARD ARM



keyboard arm
WR-3170

banana keyboard arm
WR-3170-B

- maintenance-free ball bearing track

keyboard trays/arm
color: black

25° tilt range: -15° +10°

360° swivel for easy positioning and storage

adjusts from flush to 6" below work surface

lift-n-lock height adjustment

removable memory foam palm support

*does not apply to 8491-8495 slide out tray

MONITOR ARM



monitor arm*
IN-7000-FM



monitor arm*
IN-9112-S-FM

- supplied with mounting hardware for 6 different mounting methods
- color: silver, additional colors available



desk edge
mount



through-hole
mount




wall mount



reverse wall/table
upright mount



side bolt
mount



grommet hole
mount



FLAT WHITE
color code: 248



SILVER
color code: 124



VISTA BLACK
color code: 104



PEARL WHITE
color code: 105
*only available on the 7000
monitor arm (single)

conform monitor arm color: silver

- maintenance free gas spring
- supports 5-17 lb. monitors
- integrated wire managers
- fold down for minimum desk space
- articulation range: 13"
- tilt range: -45° +90°
- clamp base

*does not apply to the 7000-fm and 9112-S-FM monitor arms

CONFORM SINGLE MONITOR ARM



monitor arm
WR-1SDA-CCG

- supports 1 computer monitor
- maximum monitor height: 16"

CONFORM DUAL MONITOR ARM



monitor arms
WR-2SDA-TPCCG

- supports 2 computer monitors
- maximum monitor height: 19.5"

TASK LIGHT



task light

- 12-5/16" wide light, 26 LEDs
- 21-5/16" wide light, 45 LEDs
- 22-9/16" wide light, 51 LEDs
- 34-3/8" wide light, 76 LEDs
- 46-3/16" wide light, 100 LEDs
- 58" wide light, 126 LEDs

- color: white
- screw mounting clips, in-line on/off rocker switch and power cord included
- 120 volts
- 4000K standard light temperature
- high efficiency with 87 Lm/W
- accessories sold separately

accessories

- 2) connection cord
- 3) direct connector
- 4) magnetic mounting clips



CPU HOLDER SUSPENDED



cpu holder

2) 8-1/4" wide - CPU29157

- color: black
- adjustable from 8-1/4"w - 3-3/4"w and 20-1/2"h - 12"h

ANTI-FATIGUE MAT



floor mat

WR-AFM-2918

- color: black
- 29.5" w x 18" d x 3/4" h
- 100% polyurethane foam with puncture resistant surface
- patent pending foot hold for convenient repositioning
- 18" beveled edges



NOVUS TOOLBAR



table size	width	surface attachment
4'	40"	clamp mount*
5'	50"	clamp mount*
6'	63"	clamp mount*

*available with cable port mount

- talented organizer for your desktop with a variety of accessories to support your workspace
- (1) aluminum slatwall element with two slots for attaching accessories, monitor arms and task lights
- (2) aluminum posts with clamp mount



NOVUS REAGENT SHELVING



table size	width	surface attachment
4'	40"	cable port mount
5'	50"	cable port mount
6'	63"	cable port mount

- talented organizer for your desktop with a variety of accessories to support your workspace
- (3) aluminum slatwall elements with two slots for attaching shelving, accessories or monitor arms
- (4) shelves; (2) shelves on the 4' option
- mount drilled through work surface



ACCESSORIES



my arm monitor arm
NV-MA-CM



LifTEC arm III
NV-LT-3



slatwall monitor arm
NV-CLII-209



multi storage board
NV-MSB-99



cell phone holder
NV-PHN-79



tablet holder
NV-TH-89



magnetic strip
NV-MS-109



pen holder
NV-PH-69



memo holder
NV-MB-39



pura-line pen cup
NV-PBH-69



penda small multi-purpose box
NV-SB-49



attenzia complete LED task light
NV-AZ-2109