

ALTUS™ TABLE SYSTEM



Responsive to Change

Altus™ tables are movable and modular – work surfaces, shelving, storage components are all independent of each other, so individual pieces can be easily and quickly reconfigured. Altus™ efficiently handles the constant change inherent in laboratories at a competitive price. Design your laboratory with adaptive and modular components and flexible layouts that can respond quickly to change.

Highly Adaptable

Flexible for changing laboratories; conveniently integrate these tables into your existing lab or use them as a basis of a freestanding easily modified new lab.

Durable - Self-supporting benches are suitable for multiple applications. They are designed to withstand heavy loads; our tables meet SEFA 10 standards.

Full Width Shelving - Allows shelves to be adjusted or removed to accommodate a wide variety of equipment.

Electrical - Pre-wired electrical outlets are conveniently located in the vertical post both above and below the worktop.

Materials, Coating Finish and Color - Altus™ table systems are available in powder coated steel, in a variety of standard or custom colors; stainless steel is an available option.

Standards - Altus™ table are available UL listed to UL 962 Standards.

Overhead Service Solutions - Altus™ tables are equipped to permit quick, easy, and affordable connection of electrical, mechanical, and data services for optimized layout flexibility. Mott offers many overhead service solutions.

Construction - Fully welded table frame construction. Single sided table is 84" high, available in multiple widths and depths.

2" x 3" Rear Posts - Allows easy access to mechanical, electrical, data and plumbing systems for servicing or change.

Plumbing - Integrated plumbing is in the vertical post, allowing for quick and easy access to services.

Data - Pre-wired data ports are available in the front post.

Cable Management - Integral cable trays are a convenient and clutter reducing solution for power and data cables.

Tailored To Your Needs - Mott has designed and manufactured thousands of tables to fit the unique requirements of every type of lab; simply contact us with your needs.



Ergonomics

Improved Productivity - The table frame can be quickly adjusted from 29" to 36" in 1" increments to easily adapt to new equipment and or users.

Height Adjustability - Raise or lower the work surface at the touch of a button with the electro-hydraulic height adjustable option. Or choose the crank-operated table frame option; legs are supplied with leveling feet.

Adjustable Shelving - Height adjustable in 1" increments with shelving options that include powder coated steel, stainless steel, wood, plastic laminate, phenolic resin or glass.



Keep It Organized

Storage Solutions - Designed to accommodate suspended, mobile and floor mounted casework; reconfigure the storage space as your tasks change. Cabinets are available in steel, stainless steel, wood, plastic laminate, phenolic resin, and steel cabinets with various options of different material fronts.

Accessories - Altus™ tables can be configured with many accessories such as keyboard trays, monitor arms, bulletin boards and shelf task lights.

Optimized Versatility

Configure Back-To-Back - Place tables back-to-back to optimize floorspace with a shared footprint; the open design is ideal for clean sightlines.

Mobile Solutions - Optional heavy-duty swivel casters have an adjustment knob for times when you need both mobility and absolute stability.

Work Table - Four-legged table with all-welded construction is height adjustable and supports heavy loads; add casters for more mobility.



Corporate Headquarters - 452 Hardy Rd, Brantford ON, Canada N3T 5L8 • T. 519.752.7825
Wood Casework Division - 562 Industrial Park Rd, WV, USA 24957 • T. 304.497.2115

www.mott.ca

QUALITY BY DESIGN

December 2020