

Acid Cabinet

Specialty Cabinet

These cabinets have been designed for the storage of acids.

Standard Features:

- Molded polyethylene liner is impervious to most acids and bases
- Radius inside corners allow for easy cleaning
- Raised lip at front edge is one inch high to contain spills
- Shelf is phenolic for chemical resistance and may be located in four different positions
- Molded liner can be removed in one piece from the front to gain access to services located behind the cabinet
- Threaded vent connection is fusion welded to the cabinet using standard plastic pipe fittings
- Door roller catch is plastic to eliminate corrosion
- Available in various heights and widths
- Interior depth of cabinet liner is 16"



mott
MANUFACTURING

ISO 9001:2008 REGISTERED

QUALITY BY DESIGN

Corporate Headquarters - 452 Hardy Rd, Brantford ON, Canada N3T 5L8 • T. 519.752.7825 • F. 519.752.2895
Wood Casework Division - 562 Industrial Park Rd, P.O. Box 156, Maxwelton, WV, USA 24957 • T. 304.497.2115

www.mott.ca