

MANCHESTER SERIES

new england lab
lab-tested furniture systems

painted steel

durable, versatile

Painted steel has long been the standard in the lab industry due to its strength and durability. Its heavy-duty steel construction is capable of supporting heavy loads. Extensive door and drawer configurations provide a variety of layout options to maximize storage in the lab. Our casework can also be constructed from satin galvanized steel for an extra level of protection against rust and corrosion in high moisture areas. The long-lasting powder coated finish, with over 20 standard color options, can also be matched to custom colors.



MANCHESTER SERIES

painted steel



Cabinet Construction

- 18 Gauge Case Construction
- Door/Drawer Fronts with Double Wall Encasing Sound Blocking Material
- Removable Backs on all Cupboard Units
- Inset Cabinet Fronts
- Welded Coved Bottoms in all Cabinets
- Optional Galvanized Steel Construction for Greater Corrosion Resistance

Cabinet Base

- 4" Integral Toe Base
- Adjustable Leveling Feet

Cabinet Shelves

- Formed Steel Shelves Adjustable on 1/2" Centers

Paint Colors

- 20 Standard Colors Available
- Two Tone Colors Available

Drawer Construction

- Type 304 Stainless Steel Drawer Body
- Painted Steel Drawer Body

Drawer Suspension

- Full Extension Ball Bearing Drawer Slides with 100 lb. Rating
- Self Closing Slides with 100 lb. Rating



Door Hardware

- Brushed Chrome 5-Knuckle Hinge
- Black 5-Knuckle Hinge
- Stainless Steel 5-Knuckle Hinge



Door/Drawer Pulls

- Aluminum Wire Pull
- Stainless Steel Wire Pull
- Brushed Aluminum "D" Pull
- Polished Chrome "D" Pull
- Black Nylon "D" Pull
- Flush Aluminum Pull
- Recessed ABS Pull
- Flush Black PVC Pull



Environmental Options

- 25.5% Post Consumer and 6.8% Pre-Consumer Recycled Steel Content