

wood veneer

natural, classic warmth

The classic look and warmth of wood has long been popular in educational and research environments. The overlay construction provides durability while still offering a pleasing appearance and has many veneer and hardware options to create the lab you desire. We offer a wide range of wood veneers from our standard oak and maple to more unique choices such as anigre, bamboo and others. Wood veneers can be stained with a variety of color options and our chemical resistant finish will withstand the rigors of your lab.



wood veneer



Cabinet Body Construction

- Veneer Core or Particle Core
- Doweled and Glued ¾" Panels
- Frameless Construction
- Cabinet Ends Finished

Cabinet Front Construction

- Full Overlay or Reveal Overlay Fronts
- ¾" Particle Core or MDF Fronts
- Built to AWI Quality Standards

Cabinet Base

- Water Resistant Plywood Base
- Veneer Core Plywood to Floor

Cabinet Shelves

- ¾" Depth with Steel Pin Shelf Supports Adjustable on 1" Centers

Drawer Construction

- Blum Metabox System
- Dovetailed Multiply Birch Box



Drawer Suspension

- Blum Metabox System
- Epoxy Coated Self-Closing, 75 lb. Slides
- Full Extension Ball Bearing, 100-150 lb. Slides

Door Hardware

- Concealed Hinge
- Brushed Chrome 5-Knuckle Hinge
- Black 5-Knuckle Hinge
- Stainless Steel 5-Knuckle Hinge
- Low Profile 3-Knuckle Hinge



Door/Drawer Pulls

- Flush Black Metal Pull
- Black Painted Wire Pull
- Aluminum Wire Pull
- Polished Chrome Recessed Pull
- Stainless Steel Wire Pull
- Brushed Aluminum Wire Pull
- Brushed Aluminum Bar Pull
- Polished Chrome Wire Pull



Edgebanding Options

- 3mm Wood Edging
- 1mm Veneer

Environmental Options

- Certified Wood Products
- No Added Urea Formaldehyde
- Recycled Content