

Delix Lab and Office Fitout - 36 Crosby Drive

OVERVIEW

- **Project:** 18,000 sq. ft. New Headquarters
- **Architect:** ACTWO Architects
- **GC:** Timberline Construction Company
- **Location:** Bedford, Massachusetts





In this new lab for Delix Lab and Office Fitout

Painted Steel Casework

Cambridge Series Benches

Epoxy Tops

Pro Series Fume Hoods

Epoxy Resin Peg Boards

Delix Therapeutics' latest Lab and Office Fit-out in Bedford, MA highlights their dedication to pushing the boundaries of neuroscience. Committed to leading the way in brain health, New England Lab's furniture and equipment were selected to ensure their new headquarters met the needs of their researchers. The lab is equipped with the adaptable Cambridge Series Benches, configured with various features like shelving, suspended cabinets, electrical raceways, and casters for flexibility and mobility. By installing fixed painted steel casework, sink stations, and fume hoods against the walls, the center of the lab allows for multiple bench layouts. For the safety of users, bench-mounted Pro Series fume hoods are provided to protect users from fumes. The lab design also includes epoxy resin pegboards, under-mount epoxy sinks, and worktops with a marine edge, enhancing the overall functionality and aesthetics of the space.