

the right formula for Alkermes- Lab Expansion

OVERVIEW

- **Project:** Laboratory Renovation
- **Architect:** Margolis + Fishman, Inc.
- **GC:** Siena Construction





In this new lab for Alkermes - Lab Expansion

Wakefield Series Designer Steel Casework

Pro Series Fume Hoods

Square - Leg Benches

Phenolic Resin Work Tops

Fully Welded "U" Frame Tables

Alkermes was looking to expand the lab of their Waltham, Massachusetts facility. Following the design details of previous projects, New England Lab provided double sided, adaptable benches with built in plumbing, electrical and data services. The flexible lab benches have painted steel suspended cabinets with maple veneer fronts, matching maple side panels and black phenolic resin work surfaces. Constant volume floor-mounted and bench-mounted fume hoods were supplied with different sash configurations. Fully welded tables and a stainless steel scullery sink with plumbing fixtures were also installed in the laboratory space.