

the right formula

for Cabot Corporation - Building 2, Labs A & B

OVERVIEW

- **Project:** Boston-based multi-laboratory research facility
- **Architect:** Margolis + Fishman, Inc.
- **GC:** Metric Construction





In this new space for Cabot Corporation - Building 2, Labs A & B

Manchester Series Painted Steel Casework

Type 304 16 Gauge Stainless Steel Tops

Harvard Series Fume Hoods

This project features two laboratories with their own color of Manchester Series painted steel casework. Both rooms required stainless steel work surfaces for the ease of cleaning. Also featured are bench and floor mounted Harvard Series fume hoods. Combination sashes with Variable Air Volume (VAV) controls keep operating costs down and provide extra protection for the user. Phenolic shelving and plumbing fixtures were also supplied.