

the right formula for Cabot Corporation - Pilot Plant

OVERVIEW

- **Project:** Large industrial facility
- **Architect:** Margolis + Fishman, Inc.
- **GC:** The Richmond Group



new england lab[®]
lab-tested furniture systems



In this new space for Cabot Corporation - Pilot Plant

Manchester Series Painted Steel Casework

Epoxy Resin Tops and Sinks

Pro Series Fume Hoods

This project features customized fume hoods utilized in a large industrial space. Several large floor-mounted hoods have a split upper combination sash and a full width lower sash. Fume hoods were also fitted with finished tops due to the highly visible space. Also featured are Manchester Series painted steel benches with epoxy resin tops and sinks.