

the right formula for Cytoc Corporation - Corporate Headquarters

OVERVIEW

- **Project:** State-of-the-art medical research facility
- **Architect:** ADD, Inc.
- **GC:** Jones Lang LaSalle (Spaulding & Slye)





In this new space for Cytoc Corporation - Corporate Headquarters

Flex Series Core Based System

Manchester Series Painted Steel Casework

Arlington Series Plastic Laminate Casework

Plastic Laminate and Epoxy Resin Tops

Harvard Series Fume Hoods

New England Lab provided laboratory furniture and equipment on a design-build basis for this project. Flex Series core system with our Manchester Series painted steel cabinets were selected to help eliminate the typical disruption of future changing needs of the company's research. Epoxy work surfaces along the length of the bench have the flexibility of being removed or height adjusted, while the cabinets are suspended or mobile for easy relocation. Specialty exhaust duct connections were built into several bench locations.