

# the right formula for Eastern Mennonite University- Suter Science Center

## OVERVIEW

- **Project:** Science Center Laboratory Renovation
- **Architect:** Blue Ridge Architects
- **GC:** Harman Construction



newenglandlab®  
lab-tested furniture systems



## In this new lab for Eastern Mennonite University- Suter Science Center

Lexington Series Wood Veneer Casework

---

RFV2 Series Fume Hoods

---

Adjustable Height Tables

---

Epoxy Resin Work Surfaces

---

Local Exhaust Arms

---

Eastern Mennonite University was looking to transform their Suter Science Center, which houses science and STEM programs, into a high tech, state-of-the-art facility. This building upgrade allows student and teacher collaboration during teaching, experimentation and other scientific activities. New England Lab provided red oak veneer casework stained with a rich dark color, as well as contrasting tan epoxy work surfaces. Teaching laboratories included group workstations with custom hydraulic, adjustable height tables adjacent to sink docking stations. Reduced face velocity fume hoods and local exhaust arms were also supplied to manage exhaust of unwanted fumes. Epoxy powder-coated galvanized steel tables were utilized in a greenhouse where moisture levels require an extra layer of corrosion resistance.