

# the right formula for Middlesex School - Centennial Science Center

## OVERVIEW

- **Project:** New and updated labs for High School
- **Architect:** Keefe Associates, Inc.
- **GC:** Linbeck





## In this new space for Middlesex School - Centennial Science Center

Lexington Series Wood Veneer Casework

---

Epoxy Resin Tops

---

Pro Series Fume Hoods

---

This high school project included chemistry, biology, physics, prep rooms, project room and control room. Lexington Series stained oak casework was used to fit out each specific lab. Special ADA height work stations and sink stations were designed within the standing high peninsula benches. Pro Series demonstration hoods allow teachers to set up and conduct experiments from within the Prep Room while students watch from the classroom. Epoxy resin tops, sinks and drying rack complete the installation.