

the right formula for Ocean Spray

OVERVIEW

- **Project:** One-of-a-kind showcase laboratory
- **Architect:** Dennis Group
- **GC:** Dennis Group



new england lab[®]
lab-tested furniture systems



In this new lab for Ocean Spray

Wakefield Series Designer Overlay Casework

Epoxy Resin Tops, Sink and Pegboards

Colortech Lab Fixtures

This project is a one of a kind showcase laboratory designed to feel more like a kitchen than a lab. To reflect Ocean Spray's classic style and building architecture, designer overlay Wakefield Series casework in bright white was fitted with white polymer door and drawer fronts that were made to look like raised panels. In the planning phase, the fronts were soaked in blueberry and cranberry concentrates to be sure they wouldn't stain. Once the raised panel fronts met the "lab-tested" standards, the lab was fully fit up with tip-out cabinets for plastic bottle recycling, under counter refrigerators with matching polymer fronts, and plenty of under-mount epoxy sinks and pegboards. Work station desk units were furnished with wood front casework and our panel based support system.