

Thermo Fisher Scientific - Patheon Greenville

OVERVIEW

- **Project:** 2nd Floor Laboratory Up-fit
- **Architect:** JacobsWyper Architects





In this new lab for Thermo Fisher Scientific - Patheon Greenville

Painted Steel Casework

Novaguard Fume Hoods

Epoxy Resin Work Tops

Safety Stations

Colortech Service Fixtures

Polypropylene Fume Hood

After the recent acquisition of Patheon, a global pharmaceutical company, Thermo Fisher Scientific wanted to upgrade their R&D, QC and QA laboratories in their facility based out of Greenville, NC. New England Lab provided the latest line of fume hoods incorporating high visibility and high performance. Alcoves with glass windows or solid panels with shelving are built into the hood side walls to provide a larger interior. Fixed islands of inset painted steel casework have epoxy resin counter tops and recessed power outlets with custom stainless steel covers. A polypropylene fume hood, powder-coated service fixtures, red acrylic splashguards and lab seating completed the space.