

Tufts University, Cummings School of Veterinary Medicine - Peabody Pavilion

OVERVIEW

- **Project:** 7,500 sq. ft. Research & Animal Care Facility
- **Architect:** Dimella Shaffer
- **Contractor:** Bowdoin Construction



new england lab
laboratory furniture systems



In this new lab for Tufts University, Cummings School of Veterinary Medicine - Peabody Pavilion

Designer Steel Overlay Casework

Optima Series Freestanding Benches

Epoxy Resin Work Tops

Cambridge Series Adjustable Height Tables

Safeguard Series Fume Hood

Tufts University recently renovated one of their buildings in North Grafton to include new science laboratories for the School of Veterinary Medicine. With updates to the facility, they will provide students with a modernized space for research. New England Lab supplied flexible freestanding bench systems that can be easily reconfigured to meet future needs. These pre-plumbed and pre-wired benches were supplied with painted steel suspended or mobile cabinets featuring plain-sliced Maple door and drawer fronts. Matching fixed cabinets were used at sink locations and where maximum storage and stability could be best utilized in the lab. Grey epoxy resin work surfaces, sinks and mobile cabinet tops were used throughout the space. Adjustable height tables, cylinder restraints, and several chemical fume hoods were also provided.