

the right formula for Wm. Jennings Bryan Dorn VAMC- Building 9 Laboratory

OVERVIEW

- **Project:** 2nd Floor Laboratory Renovation
- **Architect:** Novus Architects
- **GC:** Coburn Contractors



new england lab®
lab-tested furniture systems



In this new lab for Wm. Jennings Bryan Dorn VAMC- Building 9 Laboratory

Designer Steel Casework

Stainless Steel Casework

Plastic Laminate Casework

Observation Fume Hood

Safeguard Fume Hood

Flex Core System

The William Jennings Bryan Dorn VA Medical Center was looking to refurbish one of their research laboratories. It was essential for the lab to have flexibility so that furniture could be relocated when new equipment was required. To fit around the large columns and conceal the wires and pipes for electrical and plumbing, New England Lab supplied a 12-inch core system with cantilevered work surfaces. New England Lab featured casework made from painted steel with maple wood fronts. Stainless steel and plastic laminate casework, as well as an Observation and Safeguard fume hood were also provided to meet specific needs.