

Williams College - Science Center North Building

OVERVIEW

- **Project:** 105,000 sq. ft. New Science Center
- **Architect:** Payette
- **Contractor:** Barr & Barr
- **Location:** Williamstown, Massachusetts





In this new lab for Williams College - Science Center North Building

Designer Steel Overlay Casework

Optima Series Freestanding Benches

Epoxy Resin Work Tops

Pro Series Fume Hoods

Brightedge Series LED Task Lights

Williams College is expanding and modernizing its science facilities with the redevelopment of the Bronfman Science Center. This brand new North Building will facilitate the Geoscience, Psychology and Mathematics, and Statistics departments. New England Lab provided painted steel cabinet bodies with rift-cut White Oak fronts and gray epoxy resin work surfaces. For flexibility, freestanding benches and tables were arranged in the center of the laboratory space. Rock storage cabinets with solid wood drawers that double as portable rock trays, fill the paleo, structural geo, mineralogy, and sedimentology teaching labs. Energy-saving Brightedge LED task lights, adjustable wall shelving, and fume hoods complete the space.