

Epoxy Resin Pegboards:

All resin pegboards are manufactured out of the same material as the rest of American Epoxy Scientific's resin products and meet the same chemical and physical guidelines

- Pegboards come in a standard sizes as itemized below, However, any size pegboard can be manufactured to meet the designer's needs.
- Pegboards are fabricated with a 1/16" bevel on all exposed sides free of saw marks and have a special formulated edge dressing applied to all sides for a uniform look.
- If the back side of the pegboard is exposed it can be finished to give a smooth uniform appearance
- Double studded polypropylene pegs are supplied in 5", 6 1/2" and 8" lengths. Pegs also have a base that will help protect all glassware being stored.
- Stainless steel drip pans can also be provided to help prevent unnecessary spills in the laboratory.

WIDTH	HEIGHT	NO OF PEGS
18	24	21
20	30	34
24	24	28
24	30	38
24	36	43
30	30	51
30	36	54
36	30	60
36	36	66
42	30	68
48	30	81
48	36	88



Specifications

Thickness	1" & 3/4"
Material	Epoxy Resin
Color	Black; Charcoal Grey; Grey; Rock White; Sand; Dark Khaki; Brite White
Peg Material	Polypropylene
Peg Colors	Black, White, Grey
Peg Sizes - Lengths	5", 6-1/2", 8"
Peg Size - Diameter	0.375

