


Epoxy Resin Balance Tables:

American Epoxy Scientific LLC	1/17/2018	WEIGHT	
	MOUNTAIN HOME, AR	350 LBS.	BALANCE TABLE SET

Accessories supplied for assembly




Balance Table Bar




4 rubber pads


(8) 3 1/4" x 5/16" Hex head bolts with nuts, flat and lock washers





BALANCE TABLES CAN BE SUPPLIED IN ANY OF OUR STANDARD COLORS





All of the American Epoxy Scientific's epoxy resin balance tables are manufactured out of the same material as the rest of American Epoxy Scientific's resin products and meet the same chemical and physical guidelines .

- The balance tables are fabricated from 2-1/2" thick epoxy
- Balance tables are fabricated with a 3/16" radius on all exposed sides . The edges are also free of

saw marks and have a special formulated edge dressing applied to all sides for a uniform look.

- Standard size is 24" by 35"
- Superior load bearing capability
- Balance table bar is fabricated from steel and coated black