



DURCON EPOXY PEGBOARDS & ACCESSORIES

Durcon Epoxy Pegboards are available in 1.0" thickness in our Epoxy Signature Core colors. Pegboards are supplied with oval polypropylene pegs set at a 30° horizontal angle for effective support and drainage of lab glassware.

Optional 2" or 4" [51 mm/102mm] wide stainless steel drip trays are available for most pegboard sizes, with PVC tubing included.

PEGBOARD DIMENSIONS [IMPERIAL & METRIC]

MODEL Number	PEGS Pieces	WIDTH [inches]	HEIGHT [inches]	WIDTH [mm]	HEIGHT [mm]
1824	21	18	24	457	610
1830	25	18	30	457	762
2030	32	20	30	508	762
2424	28	24	24	610	610
2430	39	24	30	610	762
2436	36	24	36	610	914
3024	39	30	24	762	610
3030	46	30	30	762	762
3036	56	30	36	762	914
3230	53	32	30	813	762
3630	60	36	30	914	762
3636	68	36	36	914	914
4230	67	42	30	1067	762
4830	81	48	30	1219	762
4836	104	48	36	1219	914



PEG SPECIFICATIONS

MODEL Number	LENGTH [inches]
Peg Lengths	5.0", 6.5", 8.0"
Peg Diameter	0.5"
Peg Material	Polypropylene
Peg Colors	White, Black